

CLR[©]

CLINICAL LABORATORY REFERENCE

2022-2023

CLR-online.com

A SUPPLEMENT OF:



MEDICAL LABORATORY OBSERVER

FEATURING
over 4,900 tests
equipment
products
services

+

REFERENCE GUIDES
drug levels
toxicity levels
critical values
reference intervals



COLLEGE of AMERICAN
PATHOLOGISTS
Laboratory Quality Solutions

TOGETHER we move forward toward
better patient care.



Laboratory medicine keeps changing at a rapid pace. The College of American Pathologists keeps you in step with new developments—with **more than 700 proficiency testing (PT) and quality improvement programs across 16 disciplines**—so you'll have more time for what matters most: accuracy in the laboratory.

- Learn about new programs that enable you to compare multiple testing platforms for respiratory illnesses, assess your staff's technical competency for body fluid review, and more
- Get connected to direct transmission for faster and more accurate quantitative PT result reporting
- View your CAP PT results and accreditation performance with our comprehensive Performance Analytics Dashboard
- Simplify ordering with our online store—purchase PT/EQA, quality management programs, and publications

Preview our new PT programs and more at cap.org

2022-2023

CLR

CLINICAL LABORATORY REFERENCE



DEPARTMENTS

- 2** From the editor
- 4** Special Feature:
Top 10 laboratory survey findings
By Christina Wichmann

FEATURES

CLINICAL LABORATORY REFERENCE GUIDES

- 8** Table of Critical Limits (Adult)
- 9** Table of Critical Limits (Children and Newborn)
- 10** Cutoff and Toxicity Levels for Drugs-of-Abuse and Toxicology Testing
- 12** Table of Reference Intervals
- 14** Critical Values for Therapeutic Drug Levels
- 16** Pediatric Reference Table
- 38** Individual EUAs for Serology and other Adaptive Immune Response Tests for SARS-CoV-2

PRODUCT LISTINGS AND TESTS

- 18** Clinical Test and Product Information
- 26** Index of Tests, Equipment, and Services

HOW TO USE CLR

REFERENCE GUIDES

CLR is an annual supplement provided by MLO reflecting peer-reviewed clinical laboratory reference guides, as well as market resources available to clinical laboratorians.

PRODUCT INFORMATION

The product information section includes company descriptions, their essential laboratory products, and contact information for pricing and ordering.

INDEX OF TESTS, EQUIPMENT, AND SERVICES

The alphabetical index conveniently categorizes and cross-references laboratory products by test names, equipment types, and services provided.

READER FEEDBACK

Visit CLR-online.com to access the digital version. Subscribe to receive future CLR supplements here: <https://www.mlo-online.com/subscribe/>

Publishing information to help guide patient care



Medical Laboratory Observer (MLO) has been committed to producing clinical reference tables for laboratories since 1972. The 2022-23 edition of CLR celebrates 50 years of providing this valuable resource to our readers!

Critical limits and reference intervals help laboratorians interpret test results. Critical limits establish the upper and lower values that indicate when a patient's medical condition may be precarious, requiring the immediate attention of the provider who ordered the test. Reference intervals are the most common decision support tool for interpreting pathology reports, as patients' test results are compared to these ranges

to indicate normal results or those that may require clinical follow-up. These decision tools assist in one of the top goals in healthcare: timely communication.

For 15 years, I was a writer and editor of accreditation standards at The Joint Commission. In our Special Feature this month, I summarized the top 10 Joint Commission requirements and top 10 CLIA requirements that were found on laboratory surveys last year. The fifth top-scored requirement by The Joint Commission was a National Patient Safety Goal (NPSG) on improving the effectiveness of communication among healthcare staff. Specifically, this NPSG addresses the reporting of critical results of tests and diagnostic procedures on a timely basis. Having decision support tools readily accessible is a first step in preventing communication breakdowns that can contribute to missed opportunities for patient care.

As is the case every year, the editors at MLO would not have been able to produce this year's CLR issue without the help from experts such as the following:

- Gerald J. Kost, MD, PhD, MS, FACB, of the University of California-Davis Health System, for the Table of Critical Limits.
- Allison B. Chambliss, PhD, DABCC, FAACC, Director of Clinical Chemistry and Point of Care Testing, LAC and USC Medical Center and Assistant Professor of Clinical Pathology, Keck School of Medicine at USC (University of Southern California) for the Cutoff and Toxicity Levels for Drugs-of-Abuse and Toxicology Testing chart.
- S.T. Campbell, PhD, DABCC, FAACC, Department of Pathology, Montefiore Medical Center, Bronx, NY, for the Table of Reference Intervals.
- Steven W. Cotton, PhD, DABCC, FAACC, Assistant Professor in Pathology and Laboratory Medicine at the University of North Carolina, Chapel Hill, for the chart on Critical Values for Therapeutic Drug Levels.

We hope our readers find this resource beneficial. For the online version of CLR, please visit www.clr-online.com

As always, the editors at MLO are open to conversations with readers about information or suggestions to improve the charts and other content we publish in CLR. Please feel free to reach out to me at cwichmann@mlo-online.com.

Group Publisher/Executive Editor

Kristine Russell
krussell@mlo-online.com

Senior Editor

Christina Wichmann
cwichmann@mlo-online.com

Managing Editor

Gail Castanho
gcastanho@endeavorb2b.com

Graphic Artist

Patti Connors
pconnors@endeavorb2b.com

Audience Development/List Rentals

Laura Moulton
lmoulton@endeavorb2b.com

Production Director

Norma Machado
nmachado@endeavorb2b.com

ADVERTISING

East Coast/Midwest Sales (except IL)

Classified/Recruitment Advertising

Carol Vovcsko
(941) 321-2873
cvovcsko@mlo-online.com

South/West Coast/Illinois Sales

Lora Harrell
(941) 328-3707
lharrell@mlo-online.com

MLO EDITORIAL ADVISORY BOARD

John Brunstein, PhD, Biochemistry (Molecular Virology)

President & CSO
PathoID, Inc., British Columbia, Canada

John A. Gerlach, PhD, D(ABHI)

Laboratory Director
Michigan State University, East Lansing, MI

Barbara Strain, MA, SM(ASCP), CVAHP

Principal, Barbara Strain Consulting LLC
Formerly Director, Value Management
University of Virginia Health System, Charlottesville, VA

Jeffrey D. Klausner, MD, MPH

Professor of Preventive Medicine in the Division of Disease Prevention, Policy and Global Health, Department of Preventive Medicine at University of Southern California Keck School of Medicine.

Susan McQuiston, JD, MT(ASCP), SCy(ASCP)

Instructor, Biomedical Laboratory Diagnostics Program
Michigan State University, East Lansing, MI

Donna Beasley, DLM(ASCP)

Director
Huron Healthcare, Chicago, IL

Anthony Kurec, MS, H(ASCP)DLM

Clinical Associate Professor, Emeritus
SUNY Upstate Medical University, Syracuse, NY

Suzanne Butch, MLS(ASCP)CM, SBB^{CM}, DLM^{CM}

Freelance Consultant, Ann Arbor, MI

Paul R. Eden, Jr., MT(ASCP), PhD

Lt. Col., USAF (ret.)
(formerly) Chief, Laboratory Services
88th Diagnostics/Therapeutics Squadron
Wright-Patterson AFB, OH

Daniel J. Scungio, MT (ASCP), SLS, CQA (ASQ)

Consultant at Dan the Lab Safety Man and Safety Officer
at Sentara Healthcare, Norfolk, VA



CEO Chris Ferrell

CFO Mark Zadel

President June Griffin

Chief Administrative and Legal Officer Tracy Kane

COO Patrick Rains

EVP Special Projects Kristine Russell

2477 Stickney Point Rd., Suite 221B Sarasota, FL 34231

Phone: (941) 388-7050 Fax: (941) 388-7490

www.mlo-online.com



MLO - MEDICAL LABORATORY OBSERVER

(Print ISSN: 0580-7247/Online ISSN 2771-6759). Published monthly, with an additional issue in AUGUST, by Endeavor Business Media, LLC, 2477 Stickney Point Rd, Suite 221B, Sarasota, FL 34231 (941) 388-7050. Subscription rates: \$128.00/year in the U.S.; \$155.00 Canada/Mexico; Intl. subscriptions are \$221.00/year. All issues of MLO are available on microfilm from University Microfilms International, Box 78, 300 N. Zeeb Rd., Ann Arbor, MI 48106. Current single copies (if available) \$15.00 each (U.S.); and \$20.00 each (Intl.). Back issues (if available) \$1760 each (U.S.); \$22.00 each (Intl.). Payment must be made in U.S. funds on a U.S. bank/branch within the continental U.S. and accompany request. Subscription inquiries: subscriptions@endeavorb2b.com. MLO is indexed in the Cumulative Index for Nursing and Allied Health Literature and Lexis-Nexis. MLO Cover/CE, Clinical Issues, and Lab Management features are peer reviewed. Title® registered U.S. Patent Office. Copyright® 2022 by Endeavor Business Media, LLC. All rights reserved. No part of this publication may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopy, recording, or any information storage-and-retrieval system, without written permission from the publisher. Office of publication: Periodicals Postage Paid at Nashville, TN 37209 and at additional mailing offices. Postmaster: Send address changes to Omada (MLO Medical Laboratory Observer), PO Box 3257, Northbrook, IL 60065-3257. Printed in U.S.A.



DEDICATED TO IMPROVING HEALTH WORLDWIDE



QUALITY CONTROL

Leading provider of complete quality control solutions including; daily quality control, calibration verification and proficiency testing for results you can trust.



MOLECULAR DIAGNOSTICS

Molecular diagnostic solutions for clinical laboratory and POC settings. Our comprehensive test menu comprises infectious diseases including COVID-19, inherited diseases, mutation analysis and SNP genotyping.



CLINICAL CHEMISTRY ANALYZERS

Versatile range of biochemistry analyzers to suit all laboratory throughputs. Our world leading test menu delivers unrivalled consolidation and cost savings.



OPEN CHANNEL REAGENTS

The most extensive range of open channel biochemistry reagents facilitating routine and niche diagnostic testing.

Visit store.randox.com
to buy directly from Randox today



randox.com
marketing@randox.com

Product availability may vary from country to country. Some products may be for Research use Only. For more information on product application and availability, please contact your local Randox representative.



Top 10 laboratory survey findings

By Christina Wichmann

Before joining MLO, I was the editor of all the accreditation and certification standards published by The Joint Commission. I am probably one of the few people in the world who has read every single Joint Commission standard! I also helped coordinate applications and follow-up activities for deeming approval with the Centers for Medicare & Medicaid Services (CMS); in addition, I kept colleagues

apprised of regulation changes such as CLIA. This spring, both The Joint Commission and CMS published their top ten findings on laboratory surveys, which we have summarized for you below.

Table 1 summarizes the 2021 top ten findings on Joint Commission surveys. **Table 2** summarizes the top ten CLIA deficiencies found on CMS surveys. For those who may not be familiar with The Joint Commission's accreditation requirements, the standards are

an overarching statement. Under standards are a series of elements of performance, which are called EPs for short. These reflect some kind of action the organization takes and are the scorable part of a standard. CLIA regulations are written pretty broadly. More details on these regulations are on the CMS website in the State Operations Manual (SOM).¹ The SOM provides interpretive guidelines and survey procedures for each regulation. ↵

Table 1. Summary of top ten 2021 Joint Commission survey findings²

1 QSA.01.02.01, EP 2:	The laboratory conducts an investigation, performs corrective action, and maintains records of unacceptable results, late submissions, and lack of consensus in its proficiency testing program.
2 QSA.02.08.01, EP 2:	Every six months, the laboratory performs correlations to evaluate the results of the same test performed with different methodologies or instruments or at different locations.
3 HR.01.06.01, EP 18:	Staff members' competency assessments include direct observations of them performing routine patient tests, reporting test results, and performing instrument function checks and calibration.
4 QSA.02.11.01, EP 7:	As part of its quality control program, the laboratory reviews records such as equipment records and quality control summaries at least monthly.
5 NPSG.02.03.01, EP 2:	Implement procedures for reporting critical results of tests and diagnostic procedures.
6 HR.01.06.01, EP 20:	Each staff member's competence is assessed on an annual basis.
7 LD.04.05.07, EP 4:	The laboratory director or technical supervisor reviews quality control and proficiency testing data, requires corrective action on unacceptable results, and discusses issues with staff.
8 EC.02.04.03, EP 7:	The laboratory inspects, tests, and maintains all laboratory equipment.
9 HR.01.06.01, EP 3:	A qualified individual assesses a skill for staff competence.
10 QSA.05.18.01, EP 7:	The organization has policies and procedures to monitor and evaluate the patient and report suspected transfusion-related adverse events.

Table 2. Summary of top ten 2021 CLIA deficiencies from CMS surveys³

1 493.1235: Personnel Competency Assessment	As specified in the personnel requirements in subpart M, the laboratory must establish and follow written policies and procedures to assess employee and, if applicable, consultant competency.
2 493.1236(c)(1): General Lab Systems	The laboratory must define criteria for those conditions that are essential for proper storage of reagents and specimens, accurate and reliable test system operation, and test result reporting. The criteria must be consistent with the manufacturer's instructions, if provided. These conditions must be monitored and documented.
3 493.1236(c)(1): General Lab Systems	At least twice annually, the laboratory must verify the accuracy of any test or procedure it performs that is not included in subpart I or this part.
4 493.1253(b)(1): Verification of Performance Specifications	Each laboratory that introduces an unmodified, FDA-cleared or approved test system must demonstrate that it can obtain performance specifications comparable to those established by the manufacturer.
5 493.1251(b): Analytic Systems	The procedure manual must include the requirements for specimen acceptability, microscopic examination, step-by-step performance of the procedure, preparation of materials for testing, etc.
6 493.1251(a): Procedure Manual	Written procedure manual for all tests, assays, and examinations performed by the laboratory must be available to, and followed by, laboratory personnel.
7 493.1254(a)(1): Analytic Systems	Maintenance as defined by the manufacturer and with at least the frequency specified by the manufacturer.
8 493.1252(d): Analytic Systems	Reagents, solutions, culture media, control materials, calibration materials, and other supplies must not be used when they have exceeded their expiration date, have deteriorated, or are of substandard quality.
9 493.1291(c): Post Analytic Systems	The test report must indicate the following: for positive patient identification, either the patient's name and identification number, or a unique patient identifier and identification number, the name and address of the laboratory location where the test was performed, and other requirements specified in 493.1291(c).
10 493.1407(e) (4)(i): Laboratory Director's Responsibility	Ensure that the laboratory is enrolled in an HHS approved proficiency testing program for the testing performed and that the proficiency testing samples are tested as required under subpart H of this part.

REFERENCES


1. State Operations Manual. CMS.gov. Accessed June 20, 2022. https://www.cms.gov/Regulations-and-Guidance/Guidance/Manuals/downloads/som107ap_c_lab.pdf.
2. The Joint Commission. Full-year 2021 top standards noncompliance data. Joint Commission Perspectives. 2022; 42(4).
3. Centers for Medicare & Medicaid Services. Accessed June 20, 2022. <https://www.cms.gov/Regulations-and-Guidance/Legislation/CLIA/Downloads/CLIAtopten.pdf>.

Join us for the Medical Laboratory Observer's online lab management series. Each partial-day event features discussions on current topics by industry experts including videos, research, whitepapers, guides, and live sessions that keep you up-to-date with the latest information on solutions critical to the medical lab. View year-long, on-demand.




March "Technology & Trends in the Clinical Lab"

Our editors present findings from our special report on Informatics in the Laboratory, with its focus on using valuable clinical analytics to inform the decision-making process in the lab. The presentation includes additional information on valuable solutions for healthcare laboratories.



June "Lab Accreditation & MDx Workflow"

MLO presents an overview on maintaining practices that follow accreditation standards. What are the top laboratory deficiencies? How can workflow solutions help you address molecular diagnostic challenges? Do you have sufficient training in place to ensure your staff is confident and competent in their job?



November "Molecular Diagnostics for the Future Clinical Lab"

Using data gathered from our exclusive State of the Industry report on Molecular Diagnostics published in the November issue of MLO, we will explore the increasing roles of MDx and its solutions use in the clinical laboratory.

https://endeavor.swoogo.com/mlo_forum/about



ADULT

CLINICAL CHEMISTRY		LOW LIMIT		HIGH LIMIT	
Test	Units	Mean (SD)	Range	Mean (SD)	Range
Glucose	mmol/L mg/dL	2.6 (0.4) 46 (7)	1.7-3.9 30-70	26.9 (8.0) 484 (144)	6.1-55.5 110-1000
Potassium	mmol/L	2.8 (0.3)	2.5-3.6	6.2 (0.4) 8.0 (hemolyzed)	5.0-8.0
Calcium	mmol/L mg/dL	1.65 (0.17) 6.6 (0.7)	1.25-2.15 5.0-8.6	3.22 (0.22) 12.9 (0.9)	2.62-3.49 10.5-14.0
Sodium	mmol/L	120 (5)	110-137	158 (6)	145-170
CO ₂ content	mmol/L	11 (2)	5-20	40 (3)	35-50
Magnesium	mmol/L mg/dL	0.41 (0.16) 1.0 (0.4)	0.21-0.74 0.5-1.8	2.02 (0.82) 4.9 (2.0)	1.03-5.02 2.5-12.2
Phosphorus	mmol/L mg/dL	0.39 (0.10) 1.2 (0.3)	0.26-0.65 0.8-2.0	2.87 (0.48) 8.9 (1.5)	2.26-3.23 7.0-10.0
Bilirubin	μmol/L mg/dL	— —	— —	257 (86) 15 (5)	86-513 5-30
Chloride	mmol/L	75 (8)	60-90	126 (12)	115-156
Osmolality	mmol/kg	250 (13)	230-280	326 (18)	295-375
Urea nitrogen	mmol/L mg/dL	— —	— —	37.1 (21.1) 104 (59)	14.3-107.1 40-300
Uric acid	μmol/L mg/dL	— —	— —	773 (119) 13 (2)	595-892 10-15
CSF glucose	mmol/L mg/dL	2.1 (0.6) 37 (10)	1.1-2.8 20-50	24.3 (11.4) 438 (206)	13.9-38.9 250-700
Creatinine	μmol/L mg/dL	— —	— —	654 (380) 7.4 (4.3)	177-1326 2.0-15.0
Ionized calcium ⁴	mmol/L mg/dL	0.82 (0.14) 3.29 (0.56)	0.50-1.07 2.00-4.29	1.55 (0.19) 6.21 (0.76)	1.30-2.00 5.21-8.02
Lactate	mmol/L mg/dL	— —	— —	3.4 (1.3) 30.6 (11.7)	2.3-5.0 20.7-45.0

HEMATOLOGY

Hematocrit	L/L	0.18 (0.05)	0.12-0.30	0.61 (0.06)	0.54-0.80
Hemoglobin	g/L	66 (17)	40-120	199 (27)	170-300
Platelets	×10 ⁹ /L	37 (18)	10-100	910 (147)	555-1000
WBC count	×10 ⁹ /L	2.0 (0.7)	1.0-4.0	37.0 (20.7)	10.0-100.0
PT	s	—	—	27 (9)	14-40
PTT	s	—	—	68 (33)	32-150
Fibrinogen	g/L	0.88 (0.17)	0.50-1.00	7.75 (2.63)	5.00-10.00

BLOOD GASES AND PH

pCO ₂	mm Hg	19 (3)	9-25	67 (6)	50-80
pH		7.21 (0.06)	7.00-7.35	7.59 (0.03)	7.50-7.65
pO ₂	mm Hg kPa	43 (6) 5.7 (0.8)	30-55 4.0-7.3	— —	— —

Adult table modified with permission by *JAMA*, Vol. 263, pp. 704-707, 1990. CSF, cerebrospinal fluid; WBC, white blood cell; PT, prothrombin time; PTT, partial thromboplastin time. Qualitative critical results for adults¹ include the following: For *blood bank* and *immunology*—incompatible crossmatch, tests positive for syphilis (RPR or VDRL). For *microbiology* and *parasitology*—positive results from Gram stain or in culture from blood, cerebrospinal fluid, or body cavity fluid; positive India ink preparation; positive rapid antigen detection by agglutination tests for *Cryptococcus*, group B streptococci, *Haemophilus influenzae* b, or *Neisseria meningitidis*; positive results from acid-fast bacillus stain or culture; *Salmonella*, *Shigella*, or *Campylobacter* on stool culture; presence of malarial parasites. For *clinical microscopy* and *urinalysis*—elevated white blood cell count in CSF; presence of malignant cells, blasts, or microorganisms in CSF or body fluids; combination of strongly positive test results for glucose and for ketones in urine; presence of pathologic crystals (urate, cysteine, leucine, or tyrosine) on urinalysis. For *hematology*—listed frequently are the presence of blasts on blood smear; new diagnosis or findings of leukemia; presence of sickle cells (or aplastic crisis). Listed occasionally are plasma cells, band cells, atypical lymphocytes, and abnormal reticulocyte count.

Critical limits define boundaries of life-threatening values of laboratory test results. Critical results or values are those that fall outside high and low critical limits. Urgent clinician notification of critical results is the lab's responsibility. The system of critical value reporting was first implemented in a hospital by George D. Lundberg, MD, and first published in *MLO* in 1972. These tables are based on three national surveys by Gerald J. Kost, MD, PhD, MS, FACB, of the University of California-Davis Health System. Adapted with permission from his articles,^{1,4} the tables summarize critical limits used by 92 responding U.S. medical centers, including 20 trauma centers, and 39 children's hospitals. Mean and standard deviation (SD) data are presented. The frequency with which critical limits were listed can be found in the original articles.

As a rule of thumb, the "mean low" and "mean high" figures may be considered the critical limits for each test listed. Each institution should establish its own set of critical limits and clinician notification policy.

Dr. Kost conducted an independent national survey of U.S. medical centers and children's hospitals to determine ionized calcium critical limits.⁴ His extensive overview of critical limits and patient outcomes appeared in the March 1993 issue of *MLO*.³

Critical results of tests and diagnostic procedures fall significantly outside the normal range and may indicate a life-threatening situation. The objective is to provide the responsible licensed caregiver these results without delay, so the patient can be treated promptly.

The Joint Commission identifies critical values in current National Patient Safety Goals (NPSG).⁵ One goal is to report critical results of tests and diagnostic procedures on a timely basis. Inspectors check for compliance on this topic.

Elements of Performance for NPSG.02.03.01: (1) Collaborate with organization leaders to develop written procedures for managing the critical results of tests and diagnostic procedures that address the following: the definition of critical results of tests and diagnostic procedures; by whom and to whom critical results of tests and diagnostic procedures are reported; the acceptable length of time between availability and reporting of critical results of tests and diagnostic procedures; (2) implement the procedures for managing the critical results of tests and diagnostic procedures; and (3) evaluate the timeliness of reporting the critical results of tests and diagnostic procedures.

In "Global trends in critical values practices and their harmonization,"⁶ Kost and Hale investigate trends in critical values practices including improving pre-analytical processing, streamlining urgent notifications, assuring effective critical limits, assessing decision levels, and using visual logistics. Special considerations for pediatrics are addressed, since newborns/neonates must adapt to the extrauterine environment with its demands for striking physiological changes. Identifying existing personal adverse events clustered by time/location could be used to predict a patient's future adverse events. Customizing critical values is possible for some unmet needs like comparing critical values lists to national norms and clarifying protocols for repeat critical values testing. Also, site-neutral policies encourage timely

reporting, recording, and integrating critical values into a patient's closed-loop EMR.

Worldwide harmonization seems to be advancing one country at a time. Australia is moving toward harmonizing critical result management throughout the country.⁷ In Europe, the most accepted standard for accreditation and certification of clinical labs is ISO EN 15189:2012, which includes immediate notification of critical values as a special requisite. In the United States, CLSI published a new guideline.⁸ National standards of care must be considered and compared in order to harmonize critical values practices, but other than mentioning standard of care for reporting times in a table, the CLSI guideline does not adequately address, analyze, or compare standards of care in different countries.

A challenge is the harmonization of actual quantitative and qualitative triggers for emergency notifications, not just terminology. The reader can purchase GP47⁸ to learn three suggested nomenclature categories (critical-risk results, significant-risk results, and alert thresholds) and consult Appendix B therein for CAP Q-Probes critical values (renamed "alert thresholds" in a tabular summary in SI units) or access the same data free in reference⁹. However, as in recent MLO articles,¹⁰⁻¹¹ courts may not deem such Q-Probes subscriber data admissible in establishing the standard of care during litigation. Complexities of categories and how individual tests with thresholds are assigned to each of the three categories is difficult to explain to a jury.

Although controversial, repeat testing of hematology and coagulation critical values, especially in regards to pediatrics, should be noted.¹²

It is recommended that a positive test for COVID-19 be added to critical value for infectious disease.

REFERENCES

- Kost GJ. Critical limits for urgent clinician notification at U.S. medical centers. *JAMA*. 1990;263:704-707.
- Kost GJ. Critical limits for emergency clinician notification at U.S. children's hospitals. *Pediatrics*. 1991;88:597-603.
- Kost GJ. Using critical limits to improve patient outcome. *MLO*. 1993;25(3):22-27.
- Kost GJ. The significance of ionized calcium in cardiac and critical care. Availability and critical limits at U.S. medical centers and children's hospitals. *Arch Pathol Lab Med*. 1993;117:890-896.
- The Joint Commission. NPSG.02.03.01 Report critical results of tests and diagnostic procedures on a timely basis.
- Kost GJ, Hale KN. Global trends in critical values practices and their harmonization. *Clin Chem Lab Med*. 2011;49(2):167-176.
- Campbell C, Horvath A. Towards harmonization of critical laboratory result management—review of the literature and survey of Australian practices. *Clin Biochem Rev*. 2012;33(4):149-160.
- Young AB, et al. Management of critical- and significant-risk results. Wayne, PA: CLSI GP47. 2015.
- Wagar EA, et al. Critical values comparison: A College of American Pathologists Q-Probes survey of 163 clinical laboratories. *Arch Pathol Lab Med*. 2007;131:1769-1775.
- Kost GJ. Co-creating critical limits for enhanced acute care: proven need and web knowledge base. Part 1: A call to action! *MLO* 2015. Dec;47(12):34, 36-37.
- Kost GJ. Co-creating critical limits for enhanced acute care: proven need and web knowledge base. Part 2: Standard of care, what it means and how it is applied. *MLO* 2016. Jan;48(1):28-29.
- Sun SP, Garcia J, Hayden JA. Repeat critical hematology and coagulation values wastes resources, lengthens turnaround time, and delays clinical action. *Am J Clin Path*. 2018;149(3):247-252.

CHILDREN

CLINICAL CHEMISTRY		LOW LIMIT		HIGH LIMIT	
TEST	UNITS	MEAN (SD)	RANGE	MEAN (SD)	RANGE
Glucose	mmol/L	2.6 (0.5)	1.7-3.3	24.7 (8.9)	13.9-55.5
Potassium	mmol/L	2.8 (0.3)	2.0-3.5	6.4 (0.5)	5.0-8.0
Calcium	mmol/L	1.62 (0.17)	1.25-1.87	3.17 (0.22)	2.74-3.74
Sodium	mmol/L	121 (5)	110-130	156 (5)	150-170
CO ₂ content	mmol/L	11 (2)	6-18	39 (3)	33-45
Magnesium	mmol/L	0.45 (0.04)	0.41-0.49	1.77 (0.45)	1.23-3.00
Phosphorus	mmol/L	0.42 (0.16)	0.16-0.65	2.87 (0.39)	2.26-3.23
Bilirubin	μmol/L	—	—	257 (68)	86-342
Chloride	mmol/L	77 (8)	70-90	121 (5)	115-130
Osmolality	mmol/kg	253 (12)	240-270	318 (10)	300-330
Urea nitrogen	mmol/L	—	—	19.6 (11.4)	3.9-53.6
Uric acid	μmol/L	—	—	714 (119)	595-892
CSF glucose	mmol/L	1.7 (0.7)	1.1-2.8	—	—
Creatinine	μmol/L	—	—	336 (212)	221-884
Ionized calcium [†]	mmol/L	0.85 (0.13)	0.60-1.08	1.53 (0.11)	1.35-1.75
Lactate	mmol/L	—	—	4.1 (1.2)	2.4-5.5
Albumin	g/L	17 (5)	10-25	68 (10)	60-80
Ammonia	μmol/L	—	—	109 (50)	35-200
Protein	g/L	34 (5)	30-40	95 (6)	90-100
CSF protein	mg/L	—	—	1875 (854)	1000-3000

HEMATOLOGY

Hematocrit	L/L	0.20 (0.06)	0.10-0.30	0.62 (0.05)	0.54-0.70
Hemoglobin	g/L	69 (13)	50-100	208 (29)	170-250
Platelets	×10 ⁹ /L	53 (25)	20-100	916 (220)	600-1500
WBC count	×10 ⁹ /L	2.1 (0.9)	0.5-3.5	42.9 (25.1)	15.0-100.0
PT	s	—	—	21 (6)	15-35
PTT	s	—	—	62 (21)	40-100
Fibrinogen	g/L	0.77 (0.30)	0.20-12.0	—	—
Bleeding time	min	—	—	14.0 (4.0)	9.5-20.0

BLOOD GASES AND PH

pCO ₂	mm Hg	21 (6)	15-40	66 (23)	50-150
pH	—	7.21 (0.05)	7.10-7.30	7.59 (0.04)	7.50-7.70
pO ₂	mm Hg	45 (7)	30-55	124 (25)	100-150

NEWBORN

Modified potassium	CH	mmol/L	2.8 (0.4)	2.5-3.7	6.5	(See Ref. 3)
Bilirubin	CH	μmol/L	—	—	222 (86)	86-308
Hemoglobin	USMC	g/L	95 (35)	50-150	223 (23)	210-250
Hematocrit	USMC	L/L	0.33 (0.08)	0.24-0.45	0.71 (0.04)	0.65-0.75
pO ₂	USMC	mm Hg	37 (7)	30-50	92 (12)	70-100

Children and newborn tables modified with permission by *Pediatrics*, Vol. 88, pp. 597-603, 1991. CSF, cerebrospinal fluid; WBC, white blood cell; PT, prothrombin time; PTT, partial thromboplastin time; CH, Children's Hospital; USMC, U.S. Medical Centers. Qualitative critical results for children⁹ include the following: For *hematology*—presence of blasts in the blood smear; new diagnosis or findings of leukemia; presence of drepanocytes (sickle cells); atypical lymphocytes, or abnormal reticulocyte count; abnormal erythrocyte indices (mean corpuscular volume, mean corpuscular hemoglobin, mean corpuscular hemoglobin concentration). For *clinical microscopy and urinalysis*—elevated white blood cells in cCSF; presence of malignant cells, blasts, or microorganisms in CSF or body fluids; combination of strongly positive test results for glucose and for ketones in urine. For *microbiology and parasitology*—positive results from Gram stain or culture from blood, CSF, or body cavity fluid; presence of malarial parasites.

CUTOFF AND TOXICITY LEVELS FOR DRUGS-OF-ABUSE AND TOXICOLOGY TESTING

This table summarizes information for the interpretation of drugs-of-abuse toxicology assays; originally developed by the late Daniel M. Baer, MD, and updated by Richard A. Paulson, MT(ASCP), supervisor of Chemistry and Toxicology, VA Medical Center, Portland, OR. The table was updated and reviewed this year by Allison B. Chambliss, PhD, DABCC, FAACC, Director of Clinical Chemistry and Point of Care Testing, LAC and USC Medical Center, Assistant Professor of Clinical Pathology, Keck School of Medicine of USC (University of Southern California).

	Drug (and example trade names)	Common street names	Typical duration in urine after last dose	Common positive cutoff concentrations for urine screening assay*	Toxic blood level	Blood reference/therapeutic range
OPIATES	Heroin (Diacetylmorphine)	Horse, Smack, Junk, Brown Sugar, China White, H, Skag, White Horse, Skunk	1-2 days	2000 ng/mL (as morphine) 150 ng/mL (for 6-monoacetylmorphine)	>200 ng/mL	None detected
	Morphine (Duramorph)	M, White stuff, Miss Emma, Monkey	2 days	2000 ng/mL 300 ng/mL	>200 ng/mL	10-80 ng/mL
	Methadone (Dolophine)	Fizzies with MDMA Chocolate Chip Cookies, Amidone	3 days	300 ng/mL 200 ng/mL 150 ng/mL	>2000 ng/mL	For narcotic stabilization: 300-1000 ng/mL For pain: 100-400 ng/mL
	Meperidine (Demerol, Pethidine)	Demmies, Pain Killer	2-3 days	200 ng/mL	>1000 ng/mL	70-500 ng/mL
	Codeine (Analgesics with codeine)	School boy, Captain Cody, Cody, Lean, Sizzurp,	2 days	2000 ng/mL 300 ng/mL	>1000 ng/mL	10-100 ng/mL
	Tramadol ¹ (Ultram, Tramal Ultracet)	Ultra T	6 hours to 2 days	200 ng/mL	>1000 ng/mL	100-800 ng/mL, variable by patient and dosing regimen
	Oxycodone (Oxycontin, OxyIR, Percocet, Percodan)	Oxy, Oxycotton, O.C., Oxycet, Hillbilly Heroin, Percs	1-3 days	100 ng/mL 300 ng/mL	>200 ng/mL	10-100 ng/mL
	Hydrocodone (Lorcet, Vicodin, Lortab, Hycodan)	Vikes, Watson-387	1-2 days	300 ng/mL 100 ng/mL 50 ng/mL	>100 ng/mL	10-40 ng/mL
	Hydromorphone (Dilaudid)	Juice, Smack, D, Footballs, Dillies	1-2 days	2000 ng/mL 300 ng/mL	>100 ng/mL	10-30 ng/mL
	Fentanyl ² (Sublimaze, Duragesic, Actiq, Fentora)	Apache, China girl, China white, Dance fever, Friend, Goodfella, Jackpot, Murder 8, TNT, Tango and Cash	1-2 days	1 ng/mL 2 ng/mL	>3 ng/mL (naïve patients)	0.6-3.0 ng/mL (highly variable; depends on dose and route of administration)
HALLUCINOGENS	Lysergic acid diethylamide (LSD)	Acid, Blotter, Boomers, Cid, Golden Dragon, Looney Tunes, Lucy Mae, Microdots, Tabs, Yellow Sunshine	1-5 days	0.5 ng/mL 100 pg/mL	>2 ng/mL	None detected
	Marijuana and cannabinoids	Weed, Mary Jane, Ganja, Sensemilla, Blunt, Bud, Doobie, Dope, Grass, Pot, Green, Herb, Joint, Smoke, Stinkweed, Trees	Single use: 2-7 days (as Δ9-THC-COOH) Prolonged use: 1-2 months (as Δ9-THC-COOH)	15-100 ng/mL	50-200 ng/mL	None detected
	Phencyclidine	PCP, Angel dust, Hog, Embalming Fluid, Rocket Fuel, Sherms	Single use: 1 week Prolonged use: 2-4 weeks	25 ng/mL	>100 ng/mL	None detected
STIMULANTS	Cocaine	Coke, Crack, Flake, Snow	Single use: 1-3 days Prolonged use: 4 days	300 ng/mL 150 ng/mL (as metabolite benzoylecgonine)	>1000 ng/mL	None detected
	Amphetamine (Benzedrine, Dexedrine)	Speed, Bennies, Uppers, Dexies	Single use: 48 hours Prolonged use: 7-10 days	500 ng/mL 1000 ng/mL	>200 ng/mL	20-30 ng/mL
	Methylene-3,4-dioxy-Methamphetamine (MDMA)	Ecstasy, Adam, XTC, Love drug, Hug drug	Single use: 24 hours	300 ng/mL 500 ng/mL 1000 ng/mL	100-1000 ng/mL	20-30 ng/mL
	Methamphetamine (Desoxyn, Methedrine)	Speed, Meth, Crystal ice, Crank	Single use: 48 hours Prolonged use: 7-10 days	500 ng/mL 1000 ng/mL	>500 ng/mL	10-50 ng/mL

*Based on common screening assays currently in use and CAP Proficiency Testing reporting (2020) unless otherwise indicated.

Confirmation results by Gas Chromatography-Mass Spectrometry (GC-MS) or Liquid Chromatography-Mass Spectrometry/Mass Spectrometry (LC-MS/MS) vary by laboratory.

	Drug (and example trade names)	Common street names	Typical duration in urine after last dose	Common positive cutoff concentrations for urine screening assay*	Toxic blood level	Blood reference (therapeutic range)
BARBITURATES	Pentobarbital (Nembutal)	Barbs, Dolls, Phennies, Red/BlueBirds, Tooties, Yellows Yellow jackets,	2 days	300 ng/mL 200 ng/mL	>10 µg/mL	1-5 µg/mL
	Secobarbital (Seconal)	barbs, phennies, reds, red birds, yellow, yellow jacketsReds	2 days	300 ng/mL 200 ng/mL	>5 µg/mL	1-2 µg/mL
	Butabarbital (Butisol)	Goof balls, Candy, Peanuts, Stoppers	2 days	300 ng/mL 200 ng/mL	>25 µg/mL	3-25 µg/mL
	Butalbital (Fiorinal)	Goof balls, Sleepers, Stoppers, Peanuts	2 days	300 ng/mL 200 ng/mL	>20 µg/mL	5-15 µg/mL
	Phenobarbital	Barbs, phennies, reds, red birds, yellows, yellow jackets	1-3 weeks	300 ng/mL 200 ng/mL	>40 µg/mL	10-40 µg/ml
ALCOHOLS, DIOLS, & METABOLITES	Ethanol	Booze, Hooch	<1 day	10 mg/dL	80-400 mg/dL	100-150 mg/dL (for treatment of toxic alcohols)
	Methanol	Wood alcohol	<1 day	5 mg/dL (GC)	>20 mg/dL	<0.15 mg/dL
	Isopropanol	Rubbing alcohol	<1 day	5 mg/dL (GC)	>50 mg/dL	None detected
	Acetone		<1 day	5 mg/dL (GC)	>33 mg/dL	<1.0 mg/dL
	Ethylene Glycol	Antifreeze	<1 day	5 mg/dL (GC)	>50 mg/dL	None detected
SEDATIVES/HYPNOTICS/ANESTHETICS	Diazepam (Valium)	Tranks, Downers, Poles, Totem Z-bars, Zannies, Vs, Yellow/Blue Zs	Single use: Not detected Prolonged use: 5-7 days (up to 30 days)	300 ng/mL 200 ng/mL 150 ng/mL	Drug plus Metabolite: >5.0 µg/mL	Drug plus Metabolite: 0.1-1.0 µg/mL
	Oxazepam (Serax)	Tranks, Downers, Blues, Yellows,	Single use: Not detected Prolonged use: 5-7 days	300 ng/mL 200 ng/mL 150 ng/mL	>2.0 µg/mL	0.2-1.4 µg/mL
	Alprazolam (Xanax)	Tranks, Downers, Benzos, Poles, Totem Z-bars, Vs, Zannies, Yellow/Blue Zs.	Single use: Not detected Prolonged use: 5-7 days	300 ng/mL 200 ng/mL 150 ng/mL	>350 ng/mL	20-30 ng/mL
	Clonazepam (Klonopin)	Tranks, Downers, Blues, Yellows, bars, benzos, chill pills,	Single use: Not detected Prolonged use: 5-14 days	300 ng/mL 200 ng/mL 150 ng/mL	>80 ng/mL	20-70 ng/mL
	Chlordiazepoxide (Librium)	Tranks, Downers, Benzos, Poles, Totem Z-Bars, Vs, Yellow/Blue Zs, Zannies.	Single use: Not detected Prolonged use: 5-7 days	300 ng/mL 200 ng/mL 150 ng/mL	>5 µg/mL	0.7-1.0 µg/mL
	Lorazepam (Ativan, Loraz)	Tranks, Downers, Benzos, Poles, Totem Z-bars, Yellow/Blue Zs, Zannies, Vs	Single use: Not detected Prolonged use: 5-7 days	300-600 ng/mL	0.3-0.6ng/mL	50-240 ng/mL
	Flunitrazepam (Rohypnol)	Roofies, Rib, Rope, Date Rape Drug, Mexican Valium,Mind Eraser, Roaches, Roopies, Rophies	72 hours	2 ng/mL	>50 ng/mL	5-15 ng/mL
	Gamma-Hydroxybutyrate (Somatomax)	GHB, G-Caps Geebers, Fantasy, Liquid Ecstasy	12 hours	1-10 mg/L (GC; GC-MS)	>250 mg/L	48-125 mg/L (for narcolepsy)
	Ketamine Hydrochloride (Ketajet)	Special K, Lady Kay, Vitamin K, Cat Valium	<72 hours	5-10 ng/mL (GC-MS)	>7-27 µg/mL (highly variable)	0.5-5.0 µg/mL

REFERENCES

1. Morbidity and Mortality Weekly Report (MMWR). Centers for Disease Control and Prevention. <https://www.cdc.gov/mmwr/index.html>. Published June 15, 2021. Accessed June 28, 2022.
2. TIAFT Guidelines. TIAFT. <http://www.tiaft.org/tiaft-guidelines.html?setlang=en>. Accessed June 28, 2022.

TABLE OF REFERENCE INTERVALS

Specimen	Test	Conventional Units	Conversion Factor (multiply by)	SI Units
S	Albumin*	3.9-5.1 g/dL	10	39-51 g/L
B	Base excess (men)	-3.3 to +1.2 mmol/L	1	-3.3 to +1.2 mmol/L
B	Base excess (women)	-2.4 to +2.3 mmol/L	1	-2.4 to +2.3 mmol/L
P	Bicarbonate	21-29 mmol/L	1	21-29 mmol/L
S	Bilirubin, conjugated*	0.1-0.4 mg/dL	17.1	1.7-6.8 µmol/L
S	Bilirubin, total*	0.1-1.2 mg/dL	17.1	1.0-19.9 µmol/L
S/P	Calcium, total	9-10.4 mg/dL	0.25	2.24-2.6 mmol/L
B	CO ₂ content (venous)	22-26 mEq/L	1	22-26 mmol/L
S/P	Chloride*	98-107 mEq/L	1	98-107 mmol/L
S	Cholesterol (NCEP recommendation)	140-200 mg/dL	0.0259	3.6-5.2 mmol/L
S	Cortisol (a.m., total)*	5-23 µg/dL	27.6	138-635 nmol/L
S	Creatinine (Jaffe, men)*	0.9-1.3 mg/dL	88.4	80-115 µmol/L
S	Creatinine (Jaffe, women)*	0.6-1.1 mg/dL	88.4	53-97 µmol/L
S	Ferritin (men)*	39-715 ng/mL	1	39-715 µg/L
S	Ferritin (women)*	6-362 ng/mL	1	6-362 µg/L
P	Fibrinogen	200-400 mg/dL	0.01	2-4 g/L
S	Folate	9.5-39.0 ng/mL	2.265	21.5-88.4 nmol/L
S	Glucose, fasting*	74-100 mg/dL	0.0555	4.1-5.6 mmol/L
S	Haptoglobin*	30-200 mg/dL	0.01	0.3-2.0 g/L
B	Hematocrit (men)*	40.0-52.0 %	0.01	0.40-0.52 Vol fraction
B	Hematocrit (women)*	35.0-47.0 %	0.01	0.35-0.47 Vol fraction
B	Hemoglobin (men)*	14-18 g/dL	10	140-180 g/L
B	Hemoglobin (women)*	12-16 g/dL	10	120-160 g/L
S/P	Iron, total	20-168 µg/dL	0.179	3.5-30.0 µmol/L
S/P	Iron binding capacity	250-400 µg/dL	0.179	44.8-71.6 µmol/L
B	Lactate (at bed rest)	5-12 mg/dL	0.111	0.36-0.75 mmol/L
B	Lead	<25 µg/dL	0.048	<1.21 µmol/L
S	Magnesium (Atomic Absorption)	1.6-2.6 mg/dL	0.4114	0.66-1.07 mmol/L
B	MCH (RBC index)*	28.0-32.0 pg/cell	1	28.0-32.0 pg/cell
B	MCHC (RBC index)*	32.0-36.0 %	10	320-360 g/L
B	MCV (RBC index)*	83.0-95.0 fL	1	83.0-95.0 fL
S	Osmolality	280-295 mOsm/kg	1	280-295 mmol/kg
B	pCO ₂ (arterial) (men)	35-48 mm Hg	0.133	4.7-6.4 kPa
B	pCO ₂ (arterial) (women)	32-45 mm Hg	0.133	4.3-6.0 kPa
B	pH (arterial)*	7.35-7.45	1	7.35-7.45
S/P	Phosphate (as P)*	2.8-4.8 mg/dL	0.323	0.89-1.54 mmol/L
B	pO ₂ (arterial)	83-108 mm Hg	0.133	11.0-14.4 kPa
B	Platelet count	150-450 10 ³ /mm ³	1	150-450 10 ⁹ /L
S	Potassium	3.8-4.9 mEq/L	1	3.8-4.9 mmol/L
S	Protein, total (recumbent)	6.0-7.8 g/dL	10	60-78 g/L
B	RBC count (men)*	4.5-5.9 10 ⁶ /mm ³	1	4.5-5.9 10 ¹² /L
B	RBC count (women)*	4.5-5.1 10 ⁶ /mm ³	1	4.5-5.1 10 ¹² /L
S	Sodium	136-145 mEq/L	1	136-145 mmol/L
S	Thyroxine, free*	0.8-2.7 ng/dL	12.9	10.3-34.7 pmol/L
S	Thyroxine (T ₄), total (men)*	4.6-10.5 µg/dL	12.9	59-135 nmol/L
S	Thyroxine (T ₄), total (women)*	5.5-11 µg/dL	12.9	65-138 nmol/L
S	Triglyceride (NCEP recommendation)	10-150 mg/dL	0.0113	0.11-1.7 mmol/L
S	Urea nitrogen (BUN)*	8-24 mg/dL	0.357	2.8-8.6 mmol/L
S	Uric acid (men)*	3.7-7.7 mg/dL	0.059	0.22-0.46 mmol/L
S	Uric acid (women)*	2.5-6.2 mg/dL	0.059	0.15-0.37 mmol/L
S	Vitamin B12 (WHO Recommendation)	>201 pg/mL	0.733	>147 pmol/L
S	Vitamin D (25-OH)	10-65 ng/mL	2.50	25-162 nmol/L
B	WBC count	4-11 10 ³ /mm ³	1	4-11 10 ⁹ /L
S	Zinc	80-120 µg/dL	0.153	12-18 µmol/L

Specimens: B, whole blood; P, plasma; S, serum. Reference intervals depend on test method and the demographics of the normal population used.

*Adult intervals (18Y-60Y). Age specific ranges apply for pediatric and/or geriatric populations.

Source: Burtis CA, Bruns DE. *Tietz Fundamentals of Clinical Chemistry and Molecular Diagnostics*. 7th ed. St. Louis, MO; Elsevier; 2015 and Rifai, N, Horvath AR, Wittwer, CT. *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics*, 7th ed. St. Louis, MO; Elsevier; 2018 McPherson RA, Pincus MR. *Henry's Clinical Diagnosis and Management by Laboratory Methods*. 22nd ed. Philadelphia, PA: Elsevier Saunders; 22nd ed; 2011. Revised 2021 by S.T. Campbell, PhD, DABCC, FAACC, Department of Pathology, Montefiore Medical Center, Bronx, NY.



Reach for the STAR

Innovative qSTAR technology transforms molecular diagnostics and lab efficiency with fast direct amplification within 20 minutes on open RT-PCR systems.

SARS-CoV-2
RNA STAR Complete

H I G H T H R O U G H P U T / R A P I D / A C C U R A T E

This test has not been FDA cleared or approved but has been authorized by FDA for emergency use under an EUA for use by authorized laboratories; this test has been authorized only for the detection of nucleic acid from SARS-CoV-2, not for any other viruses or pathogens; and the emergency use of this test is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of *in vitro* diagnostics for detection and/or diagnosis of COVID-19 under Section 564(b)(1) of the Federal Food, Drug, and Cosmetic Act, 21 U.S.C. § 360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.

CRITICAL VALUES FOR THERAPEUTIC DRUG LEVELS

The concept of critical values for drug levels was originally developed by the late Daniel M. Baer, MD, and first published in the April 1982 issue of *MLO*. This table is an expanded version of that publication and newly revised for 2022-2023 by Steven W. Cotten PhD, DABCC, FAACC, Assistant Professor in Pathology and Laboratory Medicine, University of North Carolina at Chapel Hill.

Drug	Indication	Therapeutic Range	Critical Value	Comments
Acetaminophen	Analgesic	5-20 µg/mL	>200 µg/mL *drawn 4 hours after ingestion	*Determination if a concentration is toxic is dependent upon when it is drawn in relation to the time of ingestion of the dose. Multiple serum concentrations will be needed to monitor improvement and removal of drug.
Amikacin	Antimicrobial	Peak: 15-30 µg/mL Trough: 4-8 µg/mL	>10 µg/mL	Peak: 30 minutes after end of infusion. Trough: before next dose. Conventional dosing protocol.
Amiodarone	Antiarrhythmic	0.5-2 µg/mL	>2.5 µg/mL	Trough concentration. Serum amiodarone levels >2.5 µg/mL had a positive predictive value of 76% for adverse drug events.
Amitriptyline	Antidepressant/analgesic (neuropathic pain)	125-250 ng/mL	>500 ng/mL	Trough concentration. Life threatening cardiac toxicity and/or seizures with concentration >1000 ng/mL.
Busulfan (IV)	Anti-leukemic, Hematopoietic cell transplantation conditioning	900-1350 µmol/min	>1500 µmol/min	Area Under the Curve (AUC) calculations based on post-infusion sampling and dosing protocols vary by institution.
Carbamazepine	Antiepileptic/mood stabilizer	4-12 µg/mL	>20 µg/mL	Trough concentrations. Correlate serum concentration with clinical presentation.
Cyclosporine	Immunosuppressant	100-400 ng/mL	>500 ng/mL	Specific concentration goal dependent upon clinical situation. For concentrations drawn with intravenous therapy, blood should be drawn from site other than that where drug is infusing. (Cyclosporine adheres to plastic.) TDM levels are dependent on transplant type. Blood concentrations can be method (immunoassay or mass spectrometry) dependent.
Digoxin	Inotrope, AV node blocker	0.5-2.0 ng/mL*	>2.5 ng/mL	Samples should be drawn >8 hours after last dose. *Concentrations >1.5 ng/mL may be associated with higher mortality.
Doxepin	Antidepressant	110-250 ng/mL	>500 ng/mL	Trough concentration.
Ethosuximide	Antiepileptic	40-100 µg/mL	>200 µg/mL	Trough concentration.
Everolimus	Immunosuppressant	3-8 ng/mL	>15 ng/mL	Trough concentration. Varies by transplant protocol.
Flecainide	Antiarrhythmic	0.2-1.0 µg/mL	>1.0 µg/mL	Midpoint or trough concentration. Monitoring recommended when given concurrently with medications that may decrease metabolism (increase concentrations).
Fluconazole	Antifungal	4.0-20.0 µg/mL	None established	Limited TDM utility except in patients receiving hemodialysis.
Flucytosine	Antifungal	25-50 µg/mL	>100-200 µg/mL	Concentration should be a peak drawn 2 hours post dose.
Gentamicin	Antimicrobial	Peak: 5-10 µg/mL Trough: <2 µg/mL	Peak: >12 µg/mL Trough: >2 µg/mL	Peak: 1 hour after infusion. Trough: before next dose. Conventional dosing protocol.
Hydroxyl itraconazole	Antifungal	Not established	None established	Active metabolite of itraconazole.
Imipramine	Antidepressant	>180-240 ng/mL	>500 ng/mL	Concentration = imipramine + desipramine (metabolite).
Itraconazole	Antifungal	>0.5 µg/mL (localized) >1.0 µg/mL (systemic)	None established	Large PK variability. Should be measured within 5-7 after initiation of therapy.
Lamotrigine	Antiepileptic/mood stabilizer	1-15 µg/mL	>20 µg/mL	Trough concentration. High concentrations generally associated with increased somnolence/confusion.
Lidocaine	Antiarrhythmic	1.5-5 µg/mL	>6 µg/mL	Concentration can be drawn at any point (from separate IV line).
Lithium	Mood stabilizer	Acute: 1-1.6 mmol/L Chronic: 0.6-1.2 mmol/L	>2.0 mmol/L >5 mmol/L potentially fatal	Serum concentrations may increase in presence of hyponatremia. Concentration: 12 hours after dose.
Nortriptyline	Antidepressant/analgesic (neuropathic pain)	50-150 ng/mL	>500 ng/mL	Trough concentration.
Phenobarbital	Antiepileptic	15-40 µg/mL	>60 µg/mL	Trough concentration. Do not collect before steady state achieved.
Phenytoin	Antiepileptic	10-20 µg/mL	>20 µg/mL	Trough concentrations. Toxic >20 µg/mL (lateral nystagmus), >40 µg/mL (decreased mentation). Toxicity may occur at lower concentrations in presence of hypoalbuminemia. Consider free phenytoin.
Posaconazole	Antifungal	>0.7 µg/mL	None established	Should be measured within 7 days of initiation therapy.
Primidone	Antiepileptic	5-12 µg/mL	>15 µg/mL	Metabolized to phenobarbital.
Procainamide (PA) (metabolite: NAPA)	Antiarrhythmic	PA: 4-8 µg/mL NAPA: 10-20 µg/mL	>10 µg/mL >40 µg/mL	Mid-point or trough concentration. Procainamide monitoring is particularly important in patients who might be fast acetylators (60% to 70% of northern Europeans, and 50% of black and white Americans) and in patients with renal impairment. Procainamide and N-acetylprocainamide levels should always be measured on the same sample.
Protriptyline	Antidepressant	50-170 ng/mL	>500 ng/mL	Trough concentration.
Quinidine	Antiarrhythmic	2-5 µg/mL	>6 µg/mL	Midpoint or trough concentration.
Salicylate	Analgesic, antipyresis Anti-inflammatory	20-100 µg/mL 100-200 µg/mL	Vertigo, tinnitus 150-300 µg/mL Nausea, vomiting, hyper-ventilation 250-400 µg/mL Toxicity >500 µg/mL	Serum concentration should be used in conjunction with clinical presentation to make decision on therapy. Multiple serum concentrations will be necessary to monitor improvement and removal of drug.
Sirolimus	Immunosuppressant	4-20 ng/mL	>25 µg/mL	Trough concentration. Whole blood samples. Therapeutic levels can be lower when used in combination with other immunosuppressants. Blood concentrations can be method (immunoassay or mass spectrometry) dependent. Therapeutic levels depend on type of transplant, time post transplant, and other concomitant drug therapy.
Tacrolimus	Immunosuppressant	5-20 ng/mL	>25 ng/mL	Whole blood samples collected as trough. Therapeutic levels can be lower when used in combination with other immunosuppressants. Bias may be present between immunoassay and LC/MS methods.
Theophylline	Bronchodilator	10-20 µg/mL	>25 µg/mL	Pulmonary literature suggest that concentrations 5-15 mg/L may be as efficacious with less toxicity. Trough concentration dependent upon drug formulation.
Tobramycin	Antibacterial	Peak: 4-8 µg/mL Trough: <1.0 µg/mL	>12 µg/mL >2 µg/mL	Peak: 1 hour after end of infusion. Trough: before next dose. Conventional dosing protocol.
Valproic acid	Antiepileptic/mood stabilizer	50-125 µg/mL	>200 µg/mL	Toxicity may occur at lower concentrations in presence of hypoalbuminemia. Consider free valproic acid. Trough concentration preferred.
Vancomycin	Antimicrobial	Trough concentrations: General: 5-15 µg/mL Pneumonia: 15-20 µg/mL	Trough: >30 µg/mL	Monitoring of peaks no longer recommended. Goal trough concentration dependent upon indication. Trough: before next dose.
Voriconazole	Antifungal	1.0-5.5 µg/mL	>6 µg/mL	Should be measured within 7 days of initiation therapy.

Ranges are approximate and may vary with laboratory and/or assay. Proper interpretation of therapeutic drug concentrations requires that the specimen be drawn at an appropriate time in relation to drug administration.



C-Reactive Protein: Considered a Prognostic Biomarker to Help Manage COVID-19 Patients.¹

SEKURE® CRP Ultra Wide Range

A sensitive, ready-to-use assay for the quantitative measurement of C-reactive protein (CRP) with a wide linear range of 0.05-160 mg/L.

Measuring CRP levels provides clinical laboratories with the ability to empower healthcare providers to make more informed decisions when evaluating infections, tissue injuries, and inflammatory disorders.

The CDC Interim Guidance: "Assessment & Testing" lists CRP as a test to consider for evaluating and caring for post-COVID-19 patients.

Discover more at sekisuidiagnostics.com/crp or call 800-332-1042



© 2022 SEKISUI Diagnostics, LLC. All rights reserved. SEKURE is a registered trademark of SEKISUI Diagnostics, LLC. Because every result matters™ is a trademarks of SEKISUI Diagnostics, LLC.

¹Centers for Disease Control and Prevention (CDC). Assessment & Testing.

SEKISUI
DIAGNOSTICS
Because every result matters™

Serum Chemistries	Age	Reference Range	
Albumin*		Male (g/dL)	Female (g/dL)
	1-7 d	2.4 - 3.9	1.9 - 4.0
	8-30 d	2.1 - 4.5	1.9 - 4.4
	31-90 d	2.1 - 4.8	2.0 - 4.2
	91-180 d	2.2 - 4.9	2.3 - 4.4
	181 d - 1 y	2.2 - 4.7	2.3 - 4.7
	1 - 3 y	3.5 - 4.2	3.5 - 4.7
	4 - 6 y	3.6 - 5.2	3.6 - 5.2
	7 - 9 y	3.8 - 5.6	3.8 - 5.6
	10 - 19 y	3.8 - 5.6	3.8 - 5.6
ALT*		Male (U/L)	Female (U/L)
	1-7 d	20 - 54	21 - 54
	8-30 d	24 - 54	22 - 46
	1 - 3 mo	27 - 54	26 - 61
	4 - 6 mo	26 - 55	26 - 51
	7 - 12 mo	26 - 59	26 - 55
	1 - 3 y	19 - 59	24 - 59
	4 - 6 y	24 - 49	24 - 49
	10 - 11 y	24 - 49	24 - 44
	12 - 13 y	24 - 68	24 - 44
	14 - 15 y	24 - 59	19 - 44
	16 - 19 y	24 - 54	19 - 49
ALK*		Male (U/L)	Female (U/L)
	1-7 d	121 - 351	107 - 357
	8-30 d	138 - 486	107 - 474
	1 - 3 mo	101 - 467	125 - 547
	4 - 6 mo	94 - 425	125 - 449
	7 - 12 mo	101 - 394	101 - 431
	1 - 3 y	185 - 383	185 - 383
	4 - 6 y	191 - 450	191 - 450
	7 - 9 y	218 - 499	218 - 499
	10 - 11 y	174 - 624	169 - 657
	12 - 13 y	245 - 584	141 - 499
	14 - 15 y	169 - 618	103 - 283
	16 - 19 y	98 - 317	82 - 169
AST*		Male (U/L)	Female (U/L)
	1-7 d	26 - 98	20 - 93
	8-30 d	16 - 67	20 - 69
	1 - 3 mo	16 - 60	16 - 61
	4 - 6 mo	16 - 62	16 - 60
	7 - 12 mo	16 - 52	16 - 60
	1 - 3 y	16 - 57	16 - 57
	5 - 6 y	10 - 47	10 - 47
	7 - 9 y	10 - 36	5 - 36
	12 - 15 y	10 - 36	5 - 26
	16 - 19 y	10 - 41	0 - 26

Bilirubin, direct*	Age	Male and Female (mg/dL)	
	neonates	< 0.4	
Bilirubin, total*		Male and Female (mg/dL)	
	0 - 1 d	< 5.1	
	1 - 2 d	< 7.2	
	3 - 5 d	< 10.3	
	1 mo - adult	< 0.8	
BUN*		Male and Female (mg/dL)	
	1-7 d	1 - 13	
	8-30 d	1 - 16	
	1 - 3 mo	1 - 12	
	4 - 12 mo	1 - 14	
	1 - 3 y	4 - 17	
	4 - 13 y	6 - 17	
	14 - 19 y	7 - 21	
Calcium*		Male (mg/dL)	Female (mg/dL)
	0 - 7 d	7.6 - 11.3	7.8 - 11.2
	8 - 30 d	8.8 - 11.6	8.6 - 11.8
	31 - 90 d	8.7 - 11.2	8.2 - 11.0
	91 - 180 d	8.5 - 11.3	8.0 - 11.4
	181 - 365 d	8.0 - 10.9	8.1 - 11.0
	1 - 3 y	8.9 - 9.9	8.9 - 9.9
	4 - 11 y	9.0 - 10.1	9.0 - 10.1
	12 - 13 y	9.0 - 10.6	9.0 - 10.6
	14 - 15 y	9.3 - 10.7	9.3 - 10.7
	16 - 19 y	9.0 - 10.7	9.0 - 10.7
CO2 (venous)*		Male and Female (mmol/L)	
	0-1 wk	13 - 21	
	1 wk - 1 mo	13 - 22	
	1 - 6 mo	13 - 23	
	6 mo - 1 y	14 - 23	
	> 1 y	16 - 25	
Chloride (Cl)*		Male and Female (mmol/L)	
	0 d - 6 mo	97 - 108	
	6 mo - 1 y	97 - 106	
	> 1 y	97 - 107	
CRP*		Male (U/L)	Female (U/L)
	0 - 90 d	0.08 - 1.58	0.09 - 1.58
	91 d - 12 mo	0.08 - 1.12	0.05 - 0.79
	13 - 36 mo	0.08 - 1.12	0.08 - 0.79
	4 - 10 y	0.06 - 0.79	0.05 - 1.00
	11 - 14 y	0.08 - 0.76	0.06 - 0.81
	15 - 18 y	0.04 - 0.79	0.06 - 0.79
CK*		Male (U/L)	Female (U/L)
	0 - 90 d	29 - 303	43 - 474
	3 - 12 mo	25 - 172	27 - 242
	13 - 24 mo	28 - 162	25 - 177
	2 - 10 y	31 - 152	25 - 177
	11 - 14 y	31 - 152	31 - 172
	15 - 18 y	34 - 147	28 - 142

This chart has been updated and reviewed by Maj. Matthew Raines, MD, DABP, Medical Director and Pathologist in the U.S. Air Force.

Creatinine*		Male (mg/dL)	Female (mg/dL)
	1 - 30 d	0.5 - 1.2	0.5 - 0.9
	31 - 365 d	0.4 - 0.7	0.4 - 0.6
	1 - 3 y	0.4 - 0.7	0.4 - 0.7
	4 - 6 y	0.5 - 0.8	0.5 - 0.8
	7 - 9 y	0.6 - 0.9	0.5 - 0.9
	10 - 12 y	0.6 - 1.0	0.6 - 1.0
	13 - 15 y	0.6 - 1.2	0.7 - 1.1
	16 - 18 y	0.8 - 1.4	0.8 - 1.2
Glucose*		Male (mg/dL)	Female (mg/dL)
	0 - 1 d	36 - 110	36 - 89
	1 - 7 d	47 - 110	47 - 110
	> 7 d	54 - 117	54 - 117
Magnesium*		Male (mg/dL)	Female (mg/dL)
	0 - 90 d	1.45 - 2.15	1.49 - 2.05
	91 d - 12 mo	1.59 - 2.49	1.60 - 2.20
	13 - 36 mo	1.59 - 2.20	1.51 - 2.20
	4 - 10 y	1.49 - 2.20	1.60 - 2.50
	11 - 15 y	1.35 - 2.05	1.60 - 2.09
	16 - 18 y	1.55 - 2.10	1.49 - 1.90
Potassium (K)*		Male and Female (mmol/L)	
	0 - 1 wk	3.2 - 5.7	
	1 wk - 1 mo	3.4 - 6.2	
	1 - 6 mo	3.5 - 5.8	
	6 mo - 1 y	3.5 - 6.3	
	> 1 y	3.3 - 4.7	
Protein, total*		Male (g/dL)	Female (g/dL)
	1 - 60 d	40 - 76	3.6 - 7.0
	61 - 180 d	40 - 70	4.0 - 7.6
	181 d - 1 y	42 - 79	4.6 - 7.8
	1 - 6 y	60 - 80	6.0 - 7.8
	7 - 9 y	63 - 81	6.3 - 8.1
	10 - 19 y	64 - 86	6.4 - 8.6
Sodium (Na)*		Male and Female (mmol/L)	
	0 - 7 d	131 - 144	
	7 - 31 d	132 - 142	
	1 - 6 mo	132 - 140	
	6 mo - 1 y	131 - 140	
	> 1 y	132 - 141	
Coagulation panel PT**	Age	Male and Female (s)	
	7 - 9 y	13.1 - 15.4	
	10 - 11 y	12.9 - 15.5	
	12 - 13 y	13.1 - 15.2	
	14 - 15 y	12.9 - 15.4	
	16 - 17 y	12.6 - 15.9	
aPTT**		Male and Female (s)	
	7 - 9 y	27 - 38	
	10 - 11 y	27 - 38	
	12 - 13 y	27 - 38	
	14 - 15 y	26 - 35	
	16 - 17 y	26 - 35	

Lipid panel Cholesterol*	Age	Reference Range	
		Male (mg/dL)	Female (mg/dL)
	1 - 3 y	37 - 178	37 - 178
	4 - 6 y	103 - 184	103 - 184
	7 - 9 y	107 - 245	107 - 245
	10 - 11 y	120 - 228	122 - 242
	12 - 13 y	122 - 228	120 - 211
	14 - 15 y	101 - 222	125 - 211
	16 - 18 y	105 - 218	101 - 215
HDL-C*		Male (mg/dL)	Female (mg/dL)
	2 - < 7 y	26 - 68	16 - 62
	7 - < 12 y	28 - 76	26 - 77
	12 - < 16 y	22 - 73	28 - 79
	16 - < 19 y	28 - 72	24 - 74
LDL-C*		Male (mg/dL)	Female (mg/dL)
	13 - 36 mo	35 - 125	35 - 125
	4 - 10 y	45 - 140	35 - 135
	11 - 15 y	45 - 120	50 - 130
	16 - 18 y	55 - 120	70 - 120
Triglyceride*		Male (mg/dL)	Female (mg/dL)
	1 - 3 y	25 - 119	25 - 119
	4 - 6 y	30 - 110	30 - 110
	7 - 9 y	26 - 123	26 - 123
	10 - 11 y	22 - 131	37 - 134
	12 - 13 y	22 - 138	35 - 124
	14 - 15 y	32 - 158	36 - 129
	16 - 19 y	32 - 134	35 - 134
Iron testing Ferritin*	Age	Reference Range	
		Male (mg/dL)	Female (mg/dL)
	0 - 90 d	40 - 775	79 - 501
	91 d - 12 mo	25 - 790	25 - 560
	13 - 36 mo	12 - 501	10 - 500
	4 - 10 y	25 - 280	22 - 158
	11 - 14 y	25 - 112	15 - 112
	15 - 18 y	18 - 158	10 - 125
Iron*	Age	5 - 11 am (mcg/dL)	5 - 11 pm (mcg/dL)
	0 - 24 mo	20 - 105	20 - 140
	2 - 9 y	20 - 105	20 - 145
	10 - 14 y	20 - 100	20 - 145
	15 - 18 y	20 - 100	20 - 145
TIBC*	Age	Male (mg/dL)	Female (mg/dL)
	0 - 90 d	155 - 330	165 - 275
	91 d - 12 mo	150 - 380	250 - 455
	13 - 36 mo	215 - 420	160 - 415
	4 - 10 y	185 - 415	260 - 385
	11 - 14 y	265 - 410	250 - 420
	15 - 18 y	270 - 415	285 - 410

*Values given in this table were obtained from published studies performed on the Dimension RxL, the precursor to the current analytical platform in use at JBER Lab. Pediatric Reference Intervals, 7th ed. Washington, DC: AACC Press, 2011

**Values given in this table were obtained from published studies performed on the Stago STA-R, the precursor to the current analytical platform in use at JBER. Flanders MM, et al. Pediatric reference intervals for ten coagulation assays. Blood 2004;104:2988.

Abbreviations: ALK (alkaline phosphatase), ALT (alanine aminotransferase), aPTT (partial thromboplastin time), AST (aspartate aminotransferase), BUN (blood urea nitrogen), CK (creatinine kinase), CO2 (carbon dioxide), CRP (C-reactive protein), HDL-C (HDL cholesterol), LDL-C (LDL cholesterol), PT (prothrombin time), TIBC (total iron-binding capacity)

CLR[®]

CLINICAL LABORATORY REFERENCE

PRODUCT INFORMATION

The following section includes company descriptions with their essential laboratory products and contact information for ordering and pricing.

Also available online at CLR-online.com



FOR PROVIDERS.
BY PROVIDERS.

ACCREDITATION COMMISSION FOR HEALTH CARE (ACHC)

139 Weston Oaks Ct.
Cary, NC 27513

[Ordering/Pricing/Product Information](#)

Phone: (855) 937-2242

Website: achc.org

Email: customerservice@achc.org

A better accreditation experience

CLIA-approved since 1995, ACHC's Clinical Laboratory Accreditation program covers all specialties, including pathology, and all testing complexities. We understand labs and bring expertise and a practical, real-world focus to help solve the problems you encounter in day-to-day operations.

Transitioning is easy with ACHC's industry-leading customer service. Your dedicated Account Advisor is a process expert and will guide you at each step. Your interpretive and technical questions will be handled by a team of laboratory specialists who understand how labs work and have expertise in CLIA regulatory requirements.

Let us be your partner in achieving excellence.

See our ad on page 21



ADVANCED DATA
SYSTEMS RCM

ADVANCED DATA SYSTEMS CORPORATION

15 Prospect St
Paramus, NJ 07652

[Ordering/Pricing/Product Information](#)

Toll Free: (800) 899-4237

Phone: (800) 899-4237

Fax: (201) 497-8023

Website: www.adsc.com

Email: info@adsc.com

Ideal for claims and billing for any laboratory specialty, MedicsRCM supports a nearly 100% success rate on 1st attempt clearinghouse claims while working to ensure obtaining every dollar for every test. Payer-based verification of benefits alerts. In and out of network alerts for the No Surprises Act. Staffing issues alleviated with our +300 person outsourced workforce. Patients pay online and through their devices and call us with statement questions. Sales activity reports and an optional self serve sales portal. Financial/management reports, KPIs, analytics. LIS integrations. The same system we use, MedicsPremier from ADS, is available if in-laboratory automation is preferred on your server or our cloud.



AMERICAN PROFICIENCY INSTITUTE

1159 Business Park Dr.
Traverse City, MI 49686

[Ordering/Pricing/Product Information](#)

Phone: (800) 333-0958

Fax: (855) 900-6119

Website: www.api-pt.com

Email: customerservice@api-pt.com

American Proficiency Institute (API) offers high quality proficiency testing to over 20,000 hospital, clinic, and physician office laboratories. API's latest innovation, DataDirect, allows laboratories to upload the ability of their Laboratory Information System (LIS) to upload PT results directly to API's website. DataDirect is FREE, easy to use, and eliminates clerical errors! Join the five million results that have been successfully uploaded using this revolutionary technology. Contact us at 800-333-0958 or at DataDirect@api-pt.com and an API DataDirect Specialist will walk you through the process.

See our ad on page 29



BECKMAN COULTER, INC.

250 South Kraemer Blvd.
Brea, CA 92821

[Ordering/Pricing/Product Information](#)

Phone: (800) 526-3821

Website: www.beckmancoulter.com

Beckman Coulter is committed to advancing healthcare for every person by applying the power of science, technology and the passion and creativity of our teams to enhance the diagnostic laboratory's role in improving healthcare outcomes. Our diagnostic systems are used in complex biomedical testing, and are found in hospitals, reference laboratories and physician office settings around the globe. Beckman Coulter offers a unique combination of people, processes and solutions designed to elevate the performance of clinical laboratories and healthcare networks. We do this by accelerating care with a menu that matters, bringing the benefit of automation to all, delivering greater insights through clinical informatics and unlocking hidden value through performance partnership. An operating company of Danaher

Corporation (NYSE: DHR) since 2011, Beckman Coulter is headquartered in Brea, Calif., and has more than 11,000 global associates working diligently to make the world a healthier place.



BIO-RAD LABORATORIES

4000 Alfred Nobel Dr.

Hercules, CA 94547

[Ordering/Pricing/Product Information](#)

Toll Free: (800) 224-6723

Phone: (800) 224-6723

Fax: (510) 741-6373

Website: www.bio-rad.com/diagnostics

When it comes to quality control, Bio-Rad leads the way for clinical laboratories with a portfolio of independent quality controls and data management software solutions to facilitate effective quality testing, analysis and compliance for hundreds of the most widely used diagnostic instruments. Bio-Rad's comprehensive menu of QC products includes analytes for immunoassay, TDM, chemistry, cardiac assessment, immunology, hematology, diabetes, molecular infectious disease, serological infectious disease and more.

Additional Products & Services Provided:

- Blood Bank Testing
- Clinical Microbiology
- Diabetes and Hemoglobinopathy Testing
- HIV Testing
- Automation Platforms: BioPlex® 2200 System, EVOLIST™ microplate system.
- Infectious Disease and Autoimmune Assays
- IntelliQ
- Serum Indices



CAROLINA LIQUID CHEMISTRIES CORP.

313 Gallimore Dairy Rd
Greensboro, NC 27409

[Ordering/Pricing/Product Information](#)

Phone: (877) 722-8910

Fax: (336) 722-8915

Website: www.carolinachemistries.com

Email:

contactsales@carolinachemistries.com

Carolina Liquid Chemistries (CLC) offers high quality and low-cost chemistry analyzers and reagents. CLC offers chemistry analyzers for all sizes of laboratories and a menu of over 80

different chemistry reagents for use on a variety of instruments. The CLC floor model instruments consist of the CLC800, CLC1600, and the ultra-high speed CLC6410. CLC also offers the EasyRA benchtop analyzer that performs urine drug screens and general chemistry tests. Excellent for urgent care and oncology practices. Trade-ins offered. CLC continues to sell bar-coded reagents for use on Olympus AU instruments along with affordable service contracts. CLC is Veteran owned and ISO certified. CLC continues to offer Emergency Use Authorized (EUA) COVID-19 Rapid Antibody tests, Rapid Antigen Tests, and PCR tests.



CLINICAL SOFTWARE SOLUTIONS

PO Box 529

Queen Creek, AZ 85142

[Ordering/Pricing/Product Information](#)

Toll Free: (800) 570-0474

Phone: (480) 888-9447 Ext: 102

Fax: (480) 888-9446

Website: www.clin1mobile.net

Email: info@clin1.com

Clinical Software Solutions has provided software products and services for healthcare facilities for over 30 years. CLIN1 Laboratory and CLIN1 QC include affordable information management solutions for any type or size of lab. Intuitive screens, Auto-Release, QA, TOX Result Reviewing, Billing capture, HL7 and Web Portal options all help accomplish daily tasks with ease, increased productivity and efficiency. Unique laboratory needs? CLIN1 Web LIS has specialized hosted or on-site solutions for Clinical, Molecular, Toxicology, Cannabis testing and COVID19 testing environments available. CLIN1 Middle-ware solutions for POCT are also available. Enhance your laboratory environment with healthcare software that delivers value.

**Visit CLR online
all year at
www.clr-online.com/**



COLLEGE OF AMERICAN PATHOLOGISTS

325 Waukegan Rd
Northfield, IL 60093-2750

Ordering/Pricing/Product Information

Toll Free: (800) 323-4040

Phone: (847) 832-7000

Website: www.cap.org

Email: contactcenter@cap.org

The College of American Pathologists (CAP) is the global leader in laboratory improvement programs, including accreditation, proficiency testing/external quality assessment (PT/EQA), and continuing education programs, all of which provide laboratories with a comprehensive view of their quality processes with insight, knowledge, and peer-based educational coaching.

Supported by the scientific expertise and commitment of the world's largest board-certified pathologist organization, the CAP partners with your laboratory to foster diagnostic confidence, helping you deliver reliable test results to improve patient outcomes.

The CAP's laboratory quality solutions, initiated more than 75 years ago, provide proficiency testing to 23,000 laboratories in more than 100 countries. Drawing from the collective expertise of our pathologist members and almost 8,000 accredited laboratories worldwide, the CAP offers a comprehensive suite of programs to help you operate a high-quality laboratory and ensure that you provide accurate results with more than 700 PT/EQA programs in 16 disciplines supporting test menus with both routine and esoteric tests.

See our ad on the Inside Front Cover



COMPUGROUP MEDICAL US

10065 Red Run Blvd., Ste. 150
Owings Mills, MD 21117

Ordering/Pricing/Product Information

Phone: (800) 359-0911 Opt. 3

Website: www.cgm.com/us

Email: sales.lab.us@cgm.com

CompuGroup Medical is the leading vendor for laboratory information systems in the U.S. and offers three renowned platforms: CGM LABDAQ, CGM SCHUYLAB, and CGM AP EASY.

Designed to increase efficiency and maximize your reimbursements, CGM LABDAQ and CGM SCHUYLAB are scalable solutions that will grow with your lab. For independent labs, an optional web portal allows providers without EMR connectivity to place orders and receive results.

Similarly, CGM AP EASY simplifies the anatomic pathology workflow for AP testing labs.

CGM will guide you to the software and configuration that makes the most sense for your unique needs, workflow, and specialty.

See our ad on page 23



DIASORIN MOLECULAR

11331 Valley View Street
Cypress, CA 90630

Ordering/Pricing/Product Information

Phone: (562) 240-6500

Fax: (714) 243-4703

Website: www.molecular.diasorin.com

Email: CS.molecular@diasorin.com

DiaSorin Molecular manufactures and distributes molecular diagnostic products worldwide, helping laboratories streamline workflow and improve patient management. Our products, including our sample-to-answer Simplexa® HSV 1 & 2 and VZV kits for meningitis and encephalitis diagnostics, assist labs in the consolidation of in-house and send-out testing by simplifying processes for increased efficiency. Our LIAISON® MDX is easy-to-use and allows for convenient sample processing, thereby consolidating time and labor. Through our molecular diagnostic solutions, we seek to better patient management with fast, reliable, and clinically actionable assay results. Also visit Luminex, a DiaSorin Company, for more laboratory testing solutions.

See our ad on the Inside Back Cover

**Visit CLR online
all year at
www.clr-online.com/**



ELECTRONIC IMAGING MATERIALS

20 Forge Street
Keene, NH 03431

Ordering/Pricing/Product Information

Toll Free: (800) 535-6987

Phone: (603) 357-1459

Website: <https://barcode-labels.com/>

Email: labelexperts@eiminc.com

The Label Experts at EIM have been crafting custom, made-to-order laboratory labels since 1987. We serve both national and international clients in Laboratory, Healthcare, Pharmaceutical, Biotech, and Contract Manufacturing.

We are the creators of the premier CryoLabel® for your cryogenic needs. We specialize in durable lab labels for tubes, vials, slides, plates, microplates & cassettes, cryogenic storage, automated systems, test kits, and vaccines.

We are a Preferred Supplier of labels for Hologic® machines.

We're dubbed 'The Label Experts' for a reason and have deep level knowledge of hundreds of specialty label materials. Think it can't be labeled? Try us!



GESTALT

850 E Spokane Falls Blvd
Spokane, WA 99202

Ordering/Pricing/Product Information

Phone: (509) 492-4912

Website: www.gestaltagnostics.com

Email: info@gestaltagnostics.com

Gestalt transforms pathology through an intelligent, configurable, and vendor-neutral digital workflow that provides true interoperability enabling pathologists to diagnose diseases faster and more efficiently. Our PathFlow® solution is a cloud-based digital pathology enterprise platform that can easily be customized based on your specific preferences. This platform consists of professional, education, and research modules for ease of mixing and matching the digital needs of your facility in a single solution, freeing pathologists from tedious, repetitive, and manual

tasks allowing them to focus on their expertise, providing their invaluable experience where it matters most.



HARDY DIAGNOSTICS

1430 W McCoy Ln
Santa Maria, CA 93455

[Ordering/Pricing/Product Information](#)

Phone: (800) 266-2222

Website: hardydiagnostics.com

Email: marketing@hardydiagnostics.com

Hardy Diagnostics is an FDA-licensed manufacturer of medical devices for microbiological testing with an ISO 13485 certified Quality Management System. The company manufactures and distributes thousands of products for the culture and identification of bacteria and fungi. Hardy Diagnostics is headquartered in Santa Maria, California, and services over 10,000 laboratories across the nation. Hardy Diagnostics exports products to over 80 foreign distributors and maintains nine distribution centers nationwide. Hardy Diagnostics' mission is to produce and distribute the finest products for the detection of microorganisms and partner with its laboratory customers to diagnose and prevent disease. For more information on products and services, visit www.HardyDiagnostics.com.



HEALTH CARE LOGISTICS

P.O. Box 25
Circleville, OH 43113

[Ordering/Pricing/Product Information](#)

Toll Free: (800) 848-1633

Phone: (740) 477-3755

Fax: (800) 447-2923

Website: www.gohcl.com/

Email: hcl@GoHCL.com

Health Care Logistics® is an industry leader specializing in manufacturing, packaging and distributing unique and hard-to-find products. Top inventory lines include Lab Equipment, Biohazard Supplies and Specimen Handling Materials. The company's commitment to customer satisfaction sets it apart with services such as free samples, small package quantities and same-day shipping on most items. HCL has spent more than 40 years tailoring its products and services to the needs of its customers and will scout non-inventory items at no extra charge or manufacture custom solutions based on individual requests.

See our ad on page 40

Visit CLR online
www.clr-online.com/



IMMY

2701 Corporate Centre Drive
Norman, OK 73069

[Ordering/Pricing/Product Information](#)

Phone: (405) 360-4669

Fax: (405) 364-1058

Website: www.immy.com

Email: marketing@immy.com

IMMY focuses primarily on bringing rapid, high-quality fungal diagnostics closer to the patient. We have a long history of delivering accurate and affordable solutions that can be implemented in any laboratory setting, on any shift, around the world. IMMY is bridging the gap between fungal infections and proper treatment through rapid diagnostics for Aspergillosis, Cryptococcosis, Histoplasmosis, Coccidioidomycosis, and Blastomycosis.



KRONUS, INC.

170 South Seneca Springs Way,
Suite 105

Star, ID 83669

[Ordering/Pricing/Product Information](#)

Toll Free: (800) 457-6687

Phone: (208) 377-4800

Fax: (208) 377-4900

Website: www.kronus.com

Email: kronus@kronus.com

KRONUS continued on page 22



GET MORE FROM ACCREDITATION

ACHC is more than an accreditor. We are your partner and focus on empowering you to meet your goals. With CLIA deeming authority, we offer complete accreditation for clinical laboratories — from waived testing to pathology for inpatient, outpatient, and independent settings, including Mohs laboratories.

achc.org

(855) 937-2242 | customerservice@achc.org

PRODUCT INFORMATION

KRONUS continued from page 21

For over 30 years, KRONUS® has provided specialized immunoassay test kits to medical professionals at the world's most respected laboratory facilities. Please contact us directly with any questions, requests or comments that you may have. We welcome the opportunity to learn how KRONUS can be of assistance in meeting your laboratory's immunoassay test kit needs.

ISLET CELL AUTOIMMUNITY ZINC TRANSPORTER 8 AUTOANTIBODY

(ZnT8Ab)

- ELISA Kit, 96 Wells
- GLUTAMIC ACID DECARBOXYLASE
AUTOANTIBODY (GADAb)
- ELISA Kit, 96 Wells
- IA-2 AUTOANTIBODY (IA-2Ab)
- ELISA Kit, 96 Wells
- INSULIN AUTOANTIBODY (IAA)
- Radioimmunoassay Kit, 100 Tubes

NEUROIMMUNOLOGY

AQUAPORIN-4 AUTOANTIBODY (AQP4Ab)

- ELISA Kit, 96 Wells
- ACETYLCHOLINE RECEPTOR AUTOANTI-
BODY
- (AChRab) – BINDING; BLOCKING
- Radioimmunoassay Kits, 100 Tubes
- VOLTAGE-GATED CALCIUM CHANNEL
AUTOANTIBODY (VGCCAb): P/Q-TYPE
- Radioimmunoassay Kits, 12 and 24 Tubes

THYROID AUTOIMMUNITY

TSH RECEPTOR ANTIBODY (TRAb)

- ELISA Kit, 96 Wells
- Coated-Tube (CT) Radioimmunoassay Kits,
60 and 120 Tubes

ADRENAL AUTOIMMUNITY

STEROID 21-HYDROXYLASE ANTIBODY (21- OHAb)

- ELISA Kit, 96 Wells

FOR RESEARCH USE

3 SCREEN ISLET CELL (GAD/IA-2/ZnT8) AUTO- ANTIBODY SCREEN

- ELISA Kit, 96 Wells†

TITIN ANTIBODY

- ELISA Kit, 48 Wells†

ACETYLCHOLINE RECEPTOR ANTIBODY

(AChRab) – BINDING

- ELISA Kit, 96 Wells†

ACETYLCHOLINE RECEPTOR ANTIBODY

(AChRab) – GANGLIONIC

- RIA Kit, 25 Tubes†

VOLTAGE-GATED POTASSIUM CHANNEL

AUTOANTIBODY (VGKAb)

- Radioimmunoassay Kit, 25 Tubes†

VOLTAGE-GATED CALCIUM CHANNEL AUTO- ANTIBODY (VGCCAb): N-TYPE

- Radioimmunoassay Kit, 25 Tubes†

AQUAPORIN-4 AUTOANTIBODY (AQP4Ab)

- IFA Kit, 100 Fields†

† For Research Use Only. Not For Use in
Diagnostic Procedures.

See our ad on the Back Cover

Luminex®
A DiaSorin Company

LUMINEX

12212 Technology Blvd Ste 130
Austin, TX 78727

Ordering/Pricing/Product Information

Toll Free: (877) 785 2323
Phone: (512) 219-8020
Fax: (512) 219 5195
Website: www.luminexcorp.com
Email: orders@luminexcorp.com

At Luminex, a DiaSorin Company, our mission is to empower labs to obtain reliable, timely, and actionable answers, ultimately advancing health. We serve our customers' needs in diverse markets, including life science research, clinical diagnostics, immunology, and personalized medicine. Our goal is to globally transform healthcare by developing innovative instruments and assays that are clinically reliable, operationally efficient, and economically responsible. Check out our instruments and reagents for your targeted and syndromic molecular testing, our xMAP® Technology for your specified multiplexing needs, and our high sensitivity-imaging flow cytometry solutions. Also, visit DiaSorin Molecular for more molecular diagnostic solutions.



LUMIRADX

221 Crescent Street, 5th Floor
Waltham, MA 02453

Ordering/Pricing/Product Information

Toll Free: (888) 586-4721
Phone: (888) 586-4721
Website: www.lumiradx.com
Email: customerservices.us@lumiradx.com

LumiraDx makes diagnostic-led care simple, accessible, and affordable by providing innovative tools for both the point of care setting and large labs. Our smarter solutions are designed to enable better health outcomes for everyone. We are focused on transforming community-based healthcare by providing

critical diagnostic information to healthcare providers at the point of need, thereby enabling more informed medical decisions to improve health outcomes while lowering costs.

See our ad on page 13

meridian BIOSCIENCE®
LIFE DISCOVERED. LIFE DIAGNOSED.

MERIDIAN BIOSCIENCE INC

3471 River Hills Dr
Cincinnati, OH 45244

Ordering/Pricing/Product Information

Toll Free: (800) 543-1980
Phone: (513) 271-3700
Website: www.meridianbioscience.com
Email: mbi@meridianbioscience.com

Meridian is a fully integrated life science company that develops, manufactures, markets, and distributes a broad range of innovative diagnostic and life science products. We are dedicated to developing and delivering better solutions that give answers with speed, accuracy, and simplicity that are redefining the possibilities of life from discovery to diagnosis. Through discovery and development, we provide critical life science raw materials used in immunological and molecular tests for human, animal, plant, and environmental applications. Through diagnosis, we provide diagnostic solutions in areas including gastrointestinal and upper respiratory infections and blood lead level testing. We build relationships and provide solutions to hospitals, reference laboratories, research centers, veterinary testing centers, physician offices, diagnostics manufacturers, and biotech companies in more than 70 countries around the world.



MICROBIOLOGICS

200 Cooper Avenue North
Saint Cloud, MN 56303

Ordering/Pricing/Product Information

Toll Free: (800) 599-2847
Phone: (320) 253-7400
Website: <https://www.microbiologics.com/>
Email: info@microbiologics.com

Microbiologics is the world's leading experts and go-to collaborators in biological products and services, focused on protecting the health and safety of people around the world. For over 50 years, we have been partnering with healthcare laboratories, manufacturers and suppliers across the globe

to co-create and provide biological control materials, assay services and consulting for microbiology, molecular diagnostics and virology. With over 900 strains available in a variety of user-friendly formats, we offer the largest and most diverse line of QC microorganisms.

Targets covered include respiratory, women's health and STIs, gastrointestinal, COVID, blood culture ID, healthcare associated infections, and more.



MIRAVISTA DIAGNOSTICS

4705 Decatur Blvd.
Indianapolis, IN 46241

Ordering/Pricing/Product Information

Toll Free: (866) 647-2847

Phone: (317) 856-2681

Fax: (317) 856-3685

Website: www.miravistalabs.com

Email: info@miravistalabs.com

MiraVista Diagnostics is a specialty diagnostic laboratory and device manufacturer committed to improved patient outcomes by developing highly accurate test solutions for the diagnosis and treatment of serious fungal infections. MiraVista specializes in testing for histoplasmosis, coccidioidomycosis, blastomycosis, aspergillosis, pneumocystis and (1-3)- β -D glucan by offering antigen, antibody, and PCR testing options.

Services offered by MiraVista Diagnostics:

MVista Quantitative Antigen Detection: Histoplasma, Coccidioides, and Blastomyces EIAs

Other Antigen Detection: Aspergillus Galactomannan, B-D Glucan

Antibody Detection: MVista Histoplasma IgG IgM, MVista Coccidioides IgG IgM, Blastomyces, and Aspergillus

DNA Detection: Pneumocystis jirovecii



ORCHARD SOFTWARE CORPORATION

701 Congressional Blvd.
Carmel, IN 46032

Ordering/Pricing/Product Information

Toll Free: (800) 856-1948

Phone: (800) 856-1948

Website: www.orchardsoft.com

Email: sales@orchardsoft.com

Orchard Software Corporation is a leader in the laboratory information system industry and offers a variety of laboratory system solutions. Orchard serves more than 2,000 laboratories across the country, helping them improve efficiency, reduce errors, and enhance integration. Orchard's cloud-based solutions are installed in physician groups and clinics, hospitals, independent reference labs, student health centers, veterinary labs, public health organizations, universities, and retail facilities. Orchard provides a comprehensive suite of laboratory information systems covering clinical lab, microbiology, pathology, molecular, point-of-care testing, and toxicology.



PROLINE LABORATORY FURNITURE

10 Avco Rd.
Haverhill, MA 01835

Ordering/Pricing/Product Information

Phone: (800) 739-9067

Fax: (978) 374-4885

Website: www.1proline.com

Email: Bench@1proline.com

Proline Laboratory Furniture and tables offer a modular design that allows you to configure laboratory to meet your requirement.

Offering multiple work surface choices including but not limited to Epoxy Resin, Stainless Steel and economical Chem-Guard chemical resistant work surfaces.

Cabinets offer 90 percent extension slides with 150 pounds capacity.

Economical laboratory tables with hanging cabinets and casters offering a truly modular ergonomic design.



Puritan®

Quality since 1919

PURITAN MEDICAL PRODUCTS

31 School Street, Box 149
Guilford, ME 04443

Puritan continued on page 24



Ask about your
Buyer's Guide discount.

CGM LABDAQ

Laboratory Information System

CGM SCHUYLAB

Laboratory Information System

CGM AP EASY

Laboratory Information System



Synchronizing Healthcare



CompuGroup
Medical

Upgrade to the leading name
in laboratory software. Let
CompuGroup Medical become
your lab's trusted partner.

+1 (800) 359-0911 opt. 3
www.cgm.com/us
sales.lab.us@cgm.com

PRODUCT INFORMATION

Puritan continued from page 23

Ordering/Pricing/Product Information

Toll Free: (800) 321-2313

Phone: (207) 876-3311

Website: www.puritanmedproducts.com

Email: sales@puritanmedproducts.com

Puritan Medical Products Co., LLC, is a 103-year-old, family-owned American company known worldwide as the trusted manufacturer of single-use products for controlled environments, diagnostics and specimen collection, environmental, forensics, and genetics, medical and microbiology media industries. Puritan is FDA-registered and ISO-certified. Its products help customers collect, analyze, and transport specimens with the utmost confidence in the products' clinical integrity. Puritan manufactures all of its patented flock swabs and related products in its four US-based manufacturing facilities. For more information about Puritan Medical Products, visit www.puritanmedproducts.com

RANDOX

RANDOX LABORATORIES-US, LTD.

515 Industrial Boulevard
Kearneysville, WV 25430

Ordering/Pricing/Product Information

Toll Free: (866) 4 Randox

Phone: (304) 728-2890

Fax: (304) 728-1890

Website: www.randox.com

Email: marketing@randox.com

Randox is committed to developing revolutionary diagnostic products designed to improve health globally including chemistry, immunoassay, molecular & POC analysers in addition to high quality first and third party reagents. Leading provider of complete quality control solutions including; daily quality control, calibration verification and proficiency testing for results you can trust. Biochip Array Technology is a unique multiplexing method offering clinical, research and molecular panels. We also offer a range of COVID-19 testing solutions for all testing needs.

See our ad on page 3



Seegene TECHNOLOGIES

SEEGENE TECHNOLOGIES

325 N. Wiget Lane, Suite #140
Walnut Creek, CA 94598

Ordering/Pricing/Product Information

Phone: (925) 448-8172

Website: www.seegenetech.com

Email: info@seegenetech.com

Seegene Technologies is the US-based division of Seegene, Inc., a global developer of molecular technologies, high multiplex real-time PCR assays and Mobile laboratories.

We have one of the largest portfolios of Research Use Only (RUO) assays for infectious diseases including common respiratory infections, STD, GI, UTI, Tick-Borne Diseases and much more.

Our COVID-19 kit, Allplex™ 2019-nCoV Assay is FDA EUA. Our Novaplex™ range of assays (RUO) also include assays that are designed to detect SARS-CoV-2 variants mutations.

SEKISUI

DIAGNOSTICS

Because every result matters™

SEKISUI DIAGNOSTICS

One Wall Street
Burlington, MA 01803

Ordering/Pricing/Product Information

Phone: (781) 652-7800

Website: <https://www.sekisuidiagnostics.com>

Email: Questions@Sekisui-Dx.com

For over 40 years SEKISUI Diagnostics has been committed to providing innovative medical diagnostics to physicians and laboratories. We develop, manufacture, and supply billions of tests each year to the global healthcare market through our sales channels and our distribution/business partners.

Our product lines include:

- Clinical chemistry reagents
- Coagulation systems and reagents
- Point-of-care molecular, rapid tests and immunoassay system
- Enzymes and specialty biochemicals

See our ad on page 15

SOFTGENETICS®

Software PowerTools for Genetic Analysis

SOFTGENETICS

100 Oakwood Ave
State College, PA 16803

Ordering/Pricing/Product Information

Phone: (814) 237-9340 Ext: 202

Fax: (814) 237-9343

Website: www.softgenetics.com

Email: info@softgenetics.com

SoftGenetics LLC specializes in the development of genetic analysis software PowerTools for medical research and forensic applications. Hallmarks of SoftGenetics' tools are advanced technologies providing exceptional

accuracy, and sensitivity in an easy-to-use Windows® user interface.

Programs Include:

Mutation Surveyor® Software, for the analysis of capillary sequence data

GeneMarker® Software, for MLPA®, MSI, Fragile X, CF, Aneuploidy, all fragment analysis applications

ChimerMarker®, automated chimerism analysis and long-term monitoring software

GeneMarker®MTP, Multiple Template Processor, like having 6 GeneMarker software in one!

NextGENE® Software, analysis of NGS reads, Variants including SNVs, indels and somatic mutations, CNV, Structural Variations including gene fusions, RNA-seq

Geneticist Assistant™ Software, NGS Interpretative Workbench

GeneMarker®HTS Software, analysis of NGS mitochondrial & STR's for forensic applications

MaSTR™, Enlightened Probabilistic Mixture Analysis Software

GeneMarker®HID Software, NDIS approved Expert System for forensic STR Analysis

Visit www.softgenetics.com for technical support and to request program trials!

See our ad on page 15



WSLH Proficiency Testing

WSLH PROFICIENCY TESTING

2601 Agriculture Dr.
Madison, WI 53718

Ordering/Pricing/Product Information

Toll Free: (800) 462-5261

Phone: (608) 224-4318

Fax: (608) 265-1111

Website: <https://wslhpt.org>

Email: ptservice@slh.wisc.edu

YOUR LAB, OUR MISSION

Founded in 1966, WSLH Proficiency Testing has been providing exceptional proficiency testing services for over fifty years. WSLH Proficiency Testing (WSLH PT) is the national public health lab that is on a mission to provide reliable, affordable and easy-to-use proficiency testing products with rapid access to results. WSLH PT is approved by CMS, and accepted by CAP, COLA and The Joint Commission. Connect with our team of experienced Medical Technologists today!



LOOKING FOR A NEW CEU?



MLO-Medical Laboratory Observer has been connecting clinical Lab Directors and Managers, pathologists and oncologists with the latest product and service trends for more than 54 years. As the premier laboratory journal, our peer-reviewed content includes accredited continuing education courses with CEUs.

We have added the ability for these clinical professionals to access CEU content on this new platform. We are building our collection of courses that cover a range of topics from Molecular, Infection Prevention, Diagnostics testings, including new techniques, management and compliance. More topics to come shortly.

Our user-friendly platform makes continuing education easy and enjoyable. As you complete coursework, you'll receive a certificate as you complete your coursework that can be used to pursue your credentialing within your facilities. You'll also have a complete record of all courses you've completed.

Visit <https://educationhub.mlo-online.com/>

IN PARTNERSHIP WITH



NORTHERN ILLINOIS UNIVERSITY

College of Health and Human Sciences

PRODUCT INDEX

The following section includes company tests, equipment, products and services for the clinical lab market.

Also available online at:
CLR-online.com

#

(1,3)-Beta-D-Glucan (Fungitell[®])	
MiraVista Diagnostics.....	23
21-Hydroxylase antibody	
KRONUS, Inc.	21
21-Hydroxylase autoantibody kits	
KRONUS, Inc.	21
25-OH Vitamin D (EIA-5396) (bone health)	
Beckman Coulter, Inc.	19
6-Acetylmorphine, urine	
Carolina Liquid Chemistries Corp.....	19

A

a1-Acid glycoprotein	
Randox Laboratories-US, Ltd.	24
a1-Acid glycoprotein, immunoassay	
Randox Laboratories-US, Ltd.	24
a1-Antitrypsin assays	
Randox Laboratories-US, Ltd.	24
Accreditation services, laboratory	
Accreditation Commission for Health Care (ACHC)	18
College of American Pathologists	20
Hardy Diagnostics.....	21
Accuracy-based products	
College of American Pathologists	20
Electronic Imaging Materials	20
WSLH Proficiency Testing.....	24
Acetaminophen assays (ACETA)	
Sekisui Diagnostics.....	24
Acetaminophen, serum	
Sekisui Diagnostics.....	24
Acetylcholine receptor (AChR)-binding antibody	
KRONUS, Inc.	21
Acetylcholine receptor antibody kits	
KRONUS, Inc.	21
Acetylcholine receptor-blocking antibodies, serum	
KRONUS, Inc.	21
AchRAb (See Acetylcholine receptor antibody kits)	
KRONUS, Inc.	21
Adenovirus	
Luminex	22
Seegene Technologies.....	24
Adenovirus DNA, qualitative, RT PCR	
DiaSorin Molecular.....	20
Luminex	22
Seegene Technologies.....	24
Adenovirus testing	
Luminex	22
Seegene Technologies.....	24
Aerobic bacterial culture, general	
Hardy Diagnostics.....	21
Agglutination assays	
IMMY	21
Alanine aminotransferase assays	
Sekisui Diagnostics.....	24
Albumin	
Carolina Liquid Chemistries Corp.....	19
Randox Laboratories-US, Ltd.	24
Sekisui Diagnostics.....	24
Albumin, assays	
Randox Laboratories-US, Ltd.	24

Albumin, standard	
Randox Laboratories-US, Ltd.	24
Alcohol assays	
Carolina Liquid Chemistries Corp.....	19
Randox Laboratories-US, Ltd.	24
Aldolase	
Carolina Liquid Chemistries Corp.....	19
Randox Laboratories-US, Ltd.	24
Alkaline phosphatase (ALP)	
Carolina Liquid Chemistries Corp.....	19
Sekisui Diagnostics.....	24
Alkaline phosphatase assays	
Randox Laboratories-US, Ltd.	24
Allergen profile, latex-related	
Hardy Diagnostics.....	21
Allergy testing products and services	
Bio-Rad Laboratories	19
ALP (IFCC)	
Randox Laboratories-US, Ltd.	24
Alpha-1 acid glycoprotein (AAG)	
Randox Laboratories-US, Ltd.	24
Alpha-1-antitrypsin, serum (AAT)	
Randox Laboratories-US, Ltd.	24
Alpha-1-microglobulin	
Randox Laboratories-US, Ltd.	24
Alpha-1-microglobulin, urine	
Randox Laboratories-US, Ltd.	24
Alpha-fetoprotein	
Randox Laboratories-US, Ltd.	24
Alpha-thalassemia	
Seegene Technologies.....	24
ALT (Alanine aminotransferase)	
Carolina Liquid Chemistries Corp.....	19
Randox Laboratories-US, Ltd.	24
ALT (IFCC)	
Randox Laboratories-US, Ltd.	24
AMH (EIA-5738) (infertility)	
Beckman Coulter, Inc.	19
Amikacin assays (AMIK)	
Randox Laboratories-US, Ltd.	24
Ammonia assays (AMMON)	
Randox Laboratories-US, Ltd.	24
Ammonia, control	
Bio-Rad Laboratories	19
Randox Laboratories-US, Ltd.	24
Ammonia, plasma	
Sekisui Diagnostics.....	24
Ammonia, standard	
Randox Laboratories-US, Ltd.	24
Amphetamine assays	
Carolina Liquid Chemistries Corp.....	19
Randox Laboratories-US, Ltd.	24
Amphetamine assays, urine	
Carolina Liquid Chemistries Corp.....	19
Randox Laboratories-US, Ltd.	24
Amphetamine confirmation, urine	
Randox Laboratories-US, Ltd.	24
Amphetamine screen and confirmation, urine	
Randox Laboratories-US, Ltd.	24
Amylase (AMY)	
Carolina Liquid Chemistries Corp.....	19
Randox Laboratories-US, Ltd.	24
Amylase (AMY), serum	
Randox Laboratories-US, Ltd.	24
Amylase, assays	
Randox Laboratories-US, Ltd.	24
Sekisui Diagnostics.....	24
Anaerobic and aerobic culture	
Hardy Diagnostics.....	21

Page	Page	Page	Page
Anaerobic culture	Antisera (See Reagents, immunologic), standard	Automated systems	Randox Laboratories-US, Ltd.24
Hardy Diagnostics.....21	Randox Laboratories-US, Ltd.24	Advanced Data Systems Corporation18	Bilirubin assays, total
Anaerobic culture, aerobic culture and gram's stain	Antistreptolysin O (ASO) antibodies	Beckman Coulter, Inc.19	Randox Laboratories-US, Ltd.24
Hardy Diagnostics.....21	Randox Laboratories-US, Ltd.24	Clinical Software Solutions19	Sekisui Diagnostics.....24
Analyzer - ESR	Antistreptolysin-O assays	Hardy Diagnostics.....21	Billing and accounts receivable
Carolina Liquid Chemistries Corp.19	Randox Laboratories-US, Ltd.24	Randox Laboratories-US, Ltd.24	Advanced Data Systems Corporation18
Analyzers	Antiviral resistance, mutational analysis	Seegene Technologies.....24	Biochemistry products and services
American Proficiency Institute.....19	Seegene Technologies.....24	Randox Laboratories-US, Ltd.24	Randox Laboratories-US, Ltd.24
Beckman Coulter, Inc.19	Apo A	Autopsy	BK virus DNA, PCR
Luminex.....22	Randox Laboratories-US, Ltd.24	College of American Pathologists20	DiaSorin Molecular.....20
LumiraDx.....22	Apo B	B	Seegene Technologies.....24
Randox Laboratories-US, Ltd.24	Randox Laboratories-US, Ltd.24	Babesia microti antibody profile	Blastomyces antibodies by CF and ID
Analyzers, automated	Apolipoprotein A1 (APO A1)	DiaSorin Molecular.....20	IMMY21
American Proficiency Institute.....19	Carolina Liquid Chemistries Corp.19	Bacterial meningitis by multiplex PCR	Blastomyces antibodies, quantitative
Beckman Coulter, Inc.19	Randox Laboratories-US, Ltd.24	Seegene Technologies.....24	IMMY21
Luminex.....22	Apolipoprotein assessment	Bacterial vaginosis (sialidase), chlamydia trachomatis, neisseria gonorrhoeae, trichomonas vaginalis (NAA), vaginal yeast culture	Blastomyces antibody by immunodiffusion
Randox Laboratories-US, Ltd.24	Randox Laboratories-US, Ltd.24	Seegene Technologies.....24	IMMY21
Analyzers, benchtop	Apolipoprotein B (APO B)	Bacteriology (Microbiology)	MiraVista Diagnostics.....23
Beckman Coulter, Inc.19	Carolina Liquid Chemistries Corp.19	College of American Pathologists20	Blastomyces antigen detection
Carolina Liquid Chemistries Corp.19	Randox Laboratories-US, Ltd.24	Puritan Medical Products23	IMMY21
Hardy Diagnostics.....21	Apolipoproteins	Seegene Technologies.....24	Blastomyces antigen test
Luminex.....22	Randox Laboratories-US, Ltd.24	Barbiturate confirmation, urine	IMMY21
Randox Laboratories-US, Ltd.24	Apolipoproteins, assays	Randox Laboratories-US, Ltd.24	MiraVista Diagnostics.....23
Analyzers, blood glucose	Randox Laboratories-US, Ltd.24	Barbiturate screen, urine	Blood bank
Beckman Coulter, Inc.19	Apparel	Randox Laboratories-US, Ltd.24	Beckman Coulter, Inc.19
Analyzers, chemistry	Health Care Logistics21	Barbiturates, assays	College of American Pathologists20
Beckman Coulter, Inc.19	Arbovirus	Randox Laboratories-US, Ltd.24	Blood bank, instrumentation
Randox Laboratories-US, Ltd.24	DiaSorin Molecular.....20	Barbiturates, assays, serum	Beckman Coulter, Inc.19
Analyzers, immunoassays	Artificial Intelligence	Randox Laboratories-US, Ltd.24	Blood bank, products
Beckman Coulter, Inc.19	Gestalt.....20	Barbiturates, assays, urine	American Proficiency Institute.....19
Randox Laboratories-US, Ltd.24	ASO (See Antistreptolysin-O)	Carolina Liquid Chemistries Corp.19	Blood bank, quality control
Analyzers, semi-automated	Randox Laboratories-US, Ltd.24	Randox Laboratories-US, Ltd.24	Bio-Rad Laboratories19
Beckman Coulter, Inc.19	Aspartate aminotransferase (AST)	Barbiturates, blood	College of American Pathologists20
Randox Laboratories-US, Ltd.24	Sekisui Diagnostics.....24	Randox Laboratories-US, Ltd.24	Blood bank, reagents
Anatomic pathology	Aspergillus antigen test	Barbiturates, serum	Beckman Coulter, Inc.19
College of American Pathologists20	IMMY21	Randox Laboratories-US, Ltd.24	Bio-Rad Laboratories19
Anemia, assays	MiraVista Diagnostics.....23	Barcode products	Blood culture, routine
Beckman Coulter, Inc.19	Aspergillus antigen, EIA	Electronic Imaging Materials20	College of American Pathologists20
Antibiotic susceptibility testing (See Drug susceptibility testing)	MiraVista Diagnostics.....23	Luminex.....22	Blood gas
Hardy Diagnostics.....21	Aspergillus galactomannan, antigen detection	BCR/ABL mutation analysis	College of American Pathologists20
Antibody (See specific antibody)	IMMY21	Seegene Technologies.....24	Randox Laboratories-US, Ltd.24
Beckman Coulter, Inc.19	Aspergillus antibody by immunodiffusion	Benchmarking services	Blood glucose (See Glucose, assays)
Antibody identification	IMMY21	College of American Pathologists20	Randox Laboratories-US, Ltd.24
Beckman Coulter, Inc.19	Aspergillus antibody IgG, EIA	Benzodiazepines, assays, urine	Blood screening, reagent kits
Antibody screen	MiraVista Diagnostics.....23	Carolina Liquid Chemistries Corp.19	Randox Laboratories-US, Ltd.24
Beckman Coulter, Inc.19	Assay development service	Randox Laboratories-US, Ltd.24	Blood, culture
Antibody, purification service (See specific antibody)	Randox Laboratories-US, Ltd.24	Beta-thalassemia	Luminex.....22
Randox Laboratories-US, Ltd.24	Seegene Technologies.....24	Seegene Technologies.....24	Blood, specialty controls
Antigens	Assay matrix development service	Bicarbonate (CO2)	Randox Laboratories-US, Ltd.24
Beckman Coulter, Inc.19	Randox Laboratories-US, Ltd.24	Carolina Liquid Chemistries Corp.19	Blood, urea nitrogen (See BUN)
Hardy Diagnostics.....21	AST	Randox Laboratories-US, Ltd.24	Carolina Liquid Chemistries Corp.19
IMMY21	Carolina Liquid Chemistries Corp.19	Bile acids, assays	Bordetella parapertussis detection
Antimicrobials, assays	AST (IFCC)	Randox Laboratories-US, Ltd.24	DiaSorin Molecular.....20
Hardy Diagnostics.....21	Randox Laboratories-US, Ltd.24	Bile acids, total	Luminex.....22
Randox Laboratories-US, Ltd.24	AST (See Aspartate aminotransferase)	Carolina Liquid Chemistries Corp.19	Bordetella pertussis and B. parapertussis by PCR
Antisera (See Reagents, immunologic)	Randox Laboratories-US, Ltd.24	Randox Laboratories-US, Ltd.24	DiaSorin Molecular.....20
Randox Laboratories-US, Ltd.24	Autoimmune disease testing	Bilirubin assays, controls	Luminex.....22
Antisera (See Reagents, immunologic), bulk	KRONUS, Inc.21	Randox Laboratories-US, Ltd.24	Bordetella pertussis and bordetella parapertussis DNA, real-time PCR
Antisera (See Reagents, immunologic), polyclonal	Automated analyzers	Carolina Liquid Chemistries Corp.19	Luminex.....22
Randox Laboratories-US, Ltd.24	Beckman Coulter, Inc.19	Randox Laboratories-US, Ltd.24	Bordetella pertussis detection
Automated hematology	Hardy Diagnostics.....21	Sekisui Diagnostics.....24	DiaSorin Molecular.....20
College of American Pathologists20	Randox Laboratories-US, Ltd.24	Bilirubin assays, standard	Luminex.....22
		College of American Pathologists20	

	Page
Meridian Bioscience Inc	22
Seegene Technologies.....	24
Bordetella pertussis, nasopharyngeal culture, K	
Luminex	22
Borrelia species DNA detection by PCR (Lyme Disease)	
Seegene Technologies.....	24
BRAF V600E mutation detection by PCR	
Seegene Technologies.....	24
BUN	
Carolina Liquid Chemistries Corp.....	19
Randox Laboratories-US, Ltd.....	24
BUN (Urea)	
Carolina Liquid Chemistries Corp.....	19
BUN, standard	
Randox Laboratories-US, Ltd.....	24
Buprenorphine screen, urine	
Carolina Liquid Chemistries Corp.....	19
C	
CA-125	
Randox Laboratories-US, Ltd.....	24
CA15-3	
Randox Laboratories-US, Ltd.....	24
CA19-9	
Randox Laboratories-US, Ltd.....	24
Calcium	
Sekisui Diagnostics.....	24
Calcium assays (Electrolytes, analysis), standard	
Randox Laboratories-US, Ltd.....	24
Calcium assays (See electrolytes, analysis)	
Randox Laboratories-US, Ltd.....	24
Calcium ISE	
Randox Laboratories-US, Ltd.....	24
Calcium, serum	
Carolina Liquid Chemistries Corp.....	19
Calcium, urine	
Carolina Liquid Chemistries Corp.....	19
Calibration kits	
College of American Pathologists	20
Health Care Logistics	21
Luminex	22
Randox Laboratories-US, Ltd.....	24
WSLH Proficiency Testing.....	24
Calibration verifiers	
College of American Pathologists	20
WSLH Proficiency Testing.....	24
Calibrators	
Randox Laboratories-US, Ltd.....	24
Calibrators, custom	
Randox Laboratories-US, Ltd.....	24
Campylobacter antigen	
Meridian Bioscience Inc	22
Campylobacter detection	
Luminex	22
Seegene Technologies.....	24
Cancer (See Oncology products and services)	
Randox Laboratories-US, Ltd.....	24
Cancer (See Oncology products and services), antigens	
Randox Laboratories-US, Ltd.....	24
Cancer markers	
College of American Pathologists	20
Randox Laboratories-US, Ltd.....	24
Candida	
Luminex	22

	Page
Seegene Technologies.....	24
Candida, testing	
Luminex	22
Seegene Technologies.....	24
Cannabinoid (marijuana) confirmation, urine	
Randox Laboratories-US, Ltd.....	24
Cannabinoid (marijuana), qualitative, urine	
Carolina Liquid Chemistries Corp.....	19
Randox Laboratories-US, Ltd.....	24
Cannabinoids, assays	
Randox Laboratories-US, Ltd.....	24
Carbamazepine, serum	
Carolina Liquid Chemistries Corp.....	19
Carbon dioxide assays	
Randox Laboratories-US, Ltd.....	24
Carcinoembryonic antigen (CEA)	
Randox Laboratories-US, Ltd.....	24
Cardiac markers	
Beckman Coulter, Inc.	19
College of American Pathologists	20
Randox Laboratories-US, Ltd.....	24
Cardiac markers, controls	
Randox Laboratories-US, Ltd.....	24
Cardiac risk assessment (See Apolipoproteins; fibrinogen; Lp (a))	
Randox Laboratories-US, Ltd.....	24
Cardiac testing	
Beckman Coulter, Inc.	19
Randox Laboratories-US, Ltd.....	24
Cardiac troponin I assays	
Beckman Coulter, Inc.	19
Cardiovascular products and services	
Randox Laboratories-US, Ltd.....	24
Carts	
Health Care Logistics	21
Cell-adhesion molecule reagents	
Randox Laboratories-US, Ltd.....	24
Centrifuges	
Health Care Logistics	21
Certification	
Accreditation Commission for Health Care (ACHC)	18
College of American Pathologists	20
Chemistry	
Beckman Coulter, Inc.	19
College of American Pathologists	20
Randox Laboratories-US, Ltd.....	24
Chemistry analyzers	
Beckman Coulter, Inc.	19
Carolina Liquid Chemistries Corp.....	19
Randox Laboratories-US, Ltd.....	24
Chemistry controls	
Beckman Coulter, Inc.	19
Carolina Liquid Chemistries Corp.....	19
Randox Laboratories-US, Ltd.....	24
Chemistry products and services	
American Proficiency Institute.....	19
Beckman Coulter, Inc.	19
Carolina Liquid Chemistries Corp.....	19
College of American Pathologists	20
Randox Laboratories-US, Ltd.....	24
Chemistry test reports	
Beckman Coulter, Inc.	19
Chemistry, diagnostics	
Beckman Coulter, Inc.	19
Carolina Liquid Chemistries Corp.....	19
DiaSorin Molecular	20
Randox Laboratories-US, Ltd.....	24

	Page
Chemistry, special	
Beckman Coulter, Inc.	19
College of American Pathologists	20
Randox Laboratories-US, Ltd.....	24
Chimerism Analysis Software	
SoftGenetics.....	24
Chlamydia trachomatis by DNA probe	
DiaSorin Molecular	20
Seegene Technologies.....	24
Chlamydia trachomatis testing	
Seegene Technologies.....	24
Chloride, serum	
Carolina Liquid Chemistries Corp.....	19
Chloride, urine	
Carolina Liquid Chemistries Corp.....	19
Cholesterol (CHOL)	
Randox Laboratories-US, Ltd.....	24
Cholesterol testing	
Randox Laboratories-US, Ltd.....	24
Cholesterol, total	
Carolina Liquid Chemistries Corp.....	19
Randox Laboratories-US, Ltd.....	24
Sekisui Diagnostics.....	24
Cholinesterase (CHE)	
Carolina Liquid Chemistries Corp.....	19
CK (Creatine kinase)	
Carolina Liquid Chemistries Corp.....	19
Randox Laboratories-US, Ltd.....	24
Sekisui Diagnostics.....	24
CK (IFCC)	
Randox Laboratories-US, Ltd.....	24
CK/MB (See Creatine kinase)	
Carolina Liquid Chemistries Corp.....	19
Randox Laboratories-US, Ltd.....	24
Cleaner	
Health Care Logistics	21
Cleanrooms	
Puritan Medical Products	23
CLIA Accreditation Services	
Accreditation Commission for Health Care (ACHC)	18
Clinical chemistry, automated	
Beckman Coulter, Inc.	19
Randox Laboratories-US, Ltd.....	24
Clinical chemistry, automated, controls, reagents	
Beckman Coulter, Inc.	19
Randox Laboratories-US, Ltd.....	24
Seegene Technologies.....	24
Clinical chemistry, controls	
Beckman Coulter, Inc.	19
Bio-Rad Laboratories	19
Randox Laboratories-US, Ltd.....	24
Clinical chemistry, reagents	
Beckman Coulter, Inc.	19
Randox Laboratories-US, Ltd.....	24
Clinical microscopy	
College of American Pathologists	20
Hardy Diagnostics.....	21
Clinical pathology	
College of American Pathologists	20
Randox Laboratories-US, Ltd.....	24
Clinical toxicology	
College of American Pathologists	20
Randox Laboratories-US, Ltd.....	24
Clostridium difficile (Toxin B, cytotoxin)	
Luminex	22
Seegene Technologies.....	24
Clostridium difficile antigen (GDH)	
Meridian Bioscience Inc	22

	Page
Clostridium difficile detection (See Toxin A detection)	
Luminex	22
Seegene Technologies.....	24
Clostridium difficile detection by PCR	
Meridian Bioscience Inc	22
Seegene Technologies.....	24
Clostridium difficile toxin A assays by EIA	
Luminex	22
Clostridium difficile toxins A+B, EIA	
Meridian Bioscience Inc	22
Coagulation assays, products and services	
College of American Pathologists	20
Randox Laboratories-US, Ltd.....	24
Cocaine	
Randox Laboratories-US, Ltd.....	24
Cocaine and benzoylecgonine, blood	
Randox Laboratories-US, Ltd.....	24
Cocaine metabolite confirmation, quantitative, urine	
Randox Laboratories-US, Ltd.....	24
Cocaine metabolite screen and confirmation, blood	
Randox Laboratories-US, Ltd.....	24
Cocaine metabolite, qualitative, urine	
Carolina Liquid Chemistries Corp.....	19
Randox Laboratories-US, Ltd.....	24
Cocaine, assays	
Randox Laboratories-US, Ltd.....	24
Cocaine, metabolites, assays	
Randox Laboratories-US, Ltd.....	24
Coccidioides antigen test	
MiraVista Diagnostics.....	23
Coccidioides antibodies	
IMMY	21
Coccidioides antibodies, quantitative	
IMMY	21
MiraVista Diagnostics.....	23
Coccidioides antibody by immunodiffusion	
IMMY	21
MiraVista Diagnostics.....	23
Coccidioides immitis antibody, IgG	
IMMY	21
MiraVista Diagnostics.....	23
Coccidioides, total antibodies, ELISA	
IMMY	21
MiraVista Diagnostics.....	23
Competency (See Educational products and services)	
College of American Pathologists	20
WSLH Proficiency Testing.....	24
Competency assessment program	
Accreditation Commission for Health Care (ACHC)	18
American Proficiency Institute.....	19
College of American Pathologists	20
WSLH Proficiency Testing.....	24
Complement 3 (C3)	
Randox Laboratories-US, Ltd.....	24
Complement 4 (C4)	
Randox Laboratories-US, Ltd.....	24
Complement C2	
Randox Laboratories-US, Ltd.....	24



American
Proficiency
Institute

API DataDirect

Eliminate the #1 cause of Proficiency Testing failures,
Clerical Errors!

UPLOAD YOUR PROFICIENCY TESTING RESULTS WITHOUT AN INTERFACE!

- ▲ Automated mapping of test names
- ▲ Quantitative and Qualitative results
- ▲ No middleware or software installation
- ▲ No additional cost - it's free!

Over 5 million results have been successfully uploaded using this revolutionary technology! The process is simple:




DataDirect utilizes the ability of your LIS to run a report and create a data file which is uploaded on API's website.

Not an API customer? Take API for a Test Drive.

Sign up to receive free Chemistry PT samples and CE credits for your entire lab!

Call a Key Accounts Specialist at 800-333-0958 or email KeyAccounts@api-pt.com to enroll.

Find out more at
api-pt.com



Page	Page	Page	Page
Complement C3, serum	Randox Laboratories-US, Ltd.....24	Controls, steroids	Bio-Rad Laboratories19
Carolina Liquid Chemistries Corp.....19	Controls, autoimmune	Controls, therapeutic drug monitoring	Randox Laboratories-US, Ltd.....24
Complement C4, serum	Bio-Rad Laboratories19	Bio-Rad Laboratories19	
Carolina Liquid Chemistries Corp.....19	Controls, blood screening	Randox Laboratories-US, Ltd.....24	
Complete blood count (CBC)	Bio-Rad Laboratories19	Controls, thyronines	Bio-Rad Laboratories19
Beckman Coulter, Inc.19	Microbiologics.....22	Randox Laboratories-US, Ltd.....24	
Complete blood count (CBC) with differential (with platelet count)	Controls, cardiac marker	Controls, toxicology	Bio-Rad Laboratories19
Beckman Coulter, Inc.19	Bio-Rad Laboratories19	Randox Laboratories-US, Ltd.....24	
Complete blood count (CBC) with differential (without platelet count)	Randox Laboratories-US, Ltd.....24	Controls, tumor marker	Bio-Rad Laboratories19
Beckman Coulter, Inc.19	Controls, chemistry	Randox Laboratories-US, Ltd.....24	
Complete blood count (CBC) without differential (with platelet count)	Bio-Rad Laboratories19	Controls, urinalysis	Beckman Coulter, Inc.19
Beckman Coulter, Inc.19	Randox Laboratories-US, Ltd.....24	Bio-Rad Laboratories19	
Complete blood count (CBC) without differential (without platelet count)	Controls, coagulation assays	Randox Laboratories-US, Ltd.....24	
Beckman Coulter, Inc.19	Bio-Rad Laboratories19	Controls, urine	Beckman Coulter, Inc.19
Computational Pathology	Controls, custom	Bio-Rad Laboratories19	
Gestalt.....20	Bio-Rad Laboratories19	Randox Laboratories-US, Ltd.....24	
Computer software	Microbiologics.....22	Controls, urine chemistry	Bio-Rad Laboratories19
Advanced Data Systems	Randox Laboratories-US, Ltd.....24	Randox Laboratories-US, Ltd.....24	
Corporation18	Controls, diabetes	Controls, virology	Bio-Rad Laboratories19
Clinical Software Solutions19	Bio-Rad Laboratories19	Microbiologics22	
CompuGroup Medical US.....20	Randox Laboratories-US, Ltd.....24	Randox Laboratories-US, Ltd.....24	
Gestalt.....20	Controls, drugs of abuse	COVID-19 - Nucleic acid amplification tests (NAAT)	Luminex22
Luminex22	Bio-Rad Laboratories19	LumiraDx22	
SoftGenetics.....24	Randox Laboratories-US, Ltd.....24	Seegene Technologies.....24	
Computer systems	Controls, general chemistry	COVID-19 - Multiplex PCR Test	Carolina Liquid Chemistries Corp.....19
Advanced Data Systems	Bio-Rad Laboratories19	DiaSorin Molecular.....20	
Corporation18	Randox Laboratories-US, Ltd.....24	Luminex22	
CompuGroup Medical US.....20	Controls, hematology	Seegene Technologies.....24	
Luminex22	Beckman Coulter, Inc.19	COVID-19 - Nucleic Acid Multiplex Test	Carolina Liquid Chemistries Corp.....19
Conjugates	Bio-Rad Laboratories19	Luminex22	
Randox Laboratories-US, Ltd.....24	Randox Laboratories-US, Ltd.....24	LumiraDx22	
Conjugates, affinity purified	Controls, hormone	COVID-19 - Qualitative Assay	Beckman Coulter, Inc.19
Randox Laboratories-US, Ltd.....24	Bio-Rad Laboratories19	Carolina Liquid Chemistries Corp.....19	
Conjugates, universal	Controls, immunoassays	Luminex22	
Randox Laboratories-US, Ltd.....24	Bio-Rad Laboratories19	Seegene Technologies.....24	
Consultative and expert services	Randox Laboratories-US, Ltd.....24	COVID-19 - POC Immunoassay	Carolina Liquid Chemistries Corp.....19
Accreditation Commission for Health	Controls, immunologic	LumiraDx22	
Care (ACHC)18	Bio-Rad Laboratories19	COVID-19 - Rapid Diagnostic Test	Carolina Liquid Chemistries Corp.....19
Beckman Coulter, Inc.19	Controls, infectious disease	Hardy Diagnostics.....21	
Clinical Software Solutions19	Bio-Rad Laboratories19	LumiraDx22	
Consulting	Microbiologics.....22	Meridian Bioscience Inc22	
Accreditation Commission for Health	Randox Laboratories-US, Ltd.....24	COVID-19 - Real-time RT-PCR Test	Carolina Liquid Chemistries Corp.....19
Care (ACHC)18	Controls, ligand	DiaSorin Molecular.....20	
Beckman Coulter, Inc.19	Bio-Rad Laboratories19	Luminex22	
Clinical Software Solutions19	Controls, lipid	Seegene Technologies.....24	
Continuing education	Bio-Rad Laboratories19	COVID-19 - Serological POC Test	Carolina Liquid Chemistries Corp.....19
Accreditation Commission for Health	Controls, molecular	LumiraDx22	
Care (ACHC)18	Bio-Rad Laboratories19	C-Peptide	Beckman Coulter, Inc.19
American Proficiency Institute.....19	Microbiologics.....22	C-peptide assays	Beckman Coulter, Inc.19
College of American Pathologists20	Controls, nucleic acid	C-Reactive protein (CRP)	Carolina Liquid Chemistries Corp.....19
WSLH Proficiency Testing.....24	Bio-Rad Laboratories19	Randox Laboratories-US, Ltd.....24	
Continuing medical education	Microbiologics.....22	Sekisui Diagnostics.....24	
Accreditation Commission for Health	Controls, peptides	C-Reactive protein (CRP), quantitative	Randox Laboratories-US, Ltd.....24
Care (ACHC)18	Randox Laboratories-US, Ltd.....24	C-Reactive protein, cardiac risk (hsCRP)	Carolina Liquid Chemistries Corp.....19
College of American Pathologists20	Controls, protein	Randox Laboratories-US, Ltd.....24	
Contract manufacturing	Bio-Rad Laboratories19	C-Reactive protein, high sensitivity (CRPHS)	Randox Laboratories-US, Ltd.....24
Microbiologics.....22	Randox Laboratories-US, Ltd.....24	C-reactive protein, nephelometry kit	Randox Laboratories-US, Ltd.....24
Puritan Medical Products23	Controls, respiratory	Creatine kinase (CK) and lactic acid dehydrogenase (LD) isoenzymes	Randox Laboratories-US, Ltd.....24
Randox Laboratories-US, Ltd.....24	Bio-Rad Laboratories19	Randox Laboratories-US, Ltd.....24	
Control, quality	Microbiologics.....22	Creatine kinase (CK), MB	Randox Laboratories-US, Ltd.....24
Beckman Coulter, Inc.19	Controls, serum chemistry	Creatine kinase (CK), total, serum	Randox Laboratories-US, Ltd.....24
Bio-Rad Laboratories19	Bio-Rad Laboratories19	Creatine kinase assays (CK)	Randox Laboratories-US, Ltd.....24
Clinical Software Solutions19	Randox Laboratories-US, Ltd.....24	Creatine kinase assays, standard	Randox Laboratories-US, Ltd.....24
Microbiologics22	Controls, speciality	Creatine, serum	Randox Laboratories-US, Ltd.....24
Randox Laboratories-US, Ltd.....24	Bio-Rad Laboratories19	Creatinine (CRE)	Randox Laboratories-US, Ltd.....24
Controls	Randox Laboratories-US, Ltd.....24	Creatinine, assays	Randox Laboratories-US, Ltd.....24
Microbiologics.....22		Creatinine, standard	Randox Laboratories-US, Ltd.....24

Page	Page	Page	Page
Cytogenetics	Digital image analysis	Drug profile, urine (ten drugs) by immunoassays with GC/MS confirmation	Entamoeba detection
College of American Pathologists20	Gestalt.....20	Randox Laboratories-US, Ltd.....24	Seegene Technologies.....24
Cytokines, assays	Digital Pathology	Drug profile, urine, comprehensive (with volatiles)	Enterobacteriaceae identification tests
Randox Laboratories-US, Ltd.....24	Gestalt.....20	Randox Laboratories-US, Ltd.....24	Hardy Diagnostics.....21
Cytology products and services	Digitoxin, assays and therapeutic drug monitoring	Drug profile, urine, comprehensive (without volatiles)	Enterohemorrhagic E. coli (EHEC) culture
College of American Pathologists20	Randox Laboratories-US, Ltd.....24	Randox Laboratories-US, Ltd.....24	Hardy Diagnostics.....21
Electronic Imaging Materials20	Digitoxin, serum	Drug susceptibility testing	Enterohemorrhagic E. coli (EHEC) shiga toxin, EIA
Puritan Medical Products23	Randox Laboratories-US, Ltd.....24	College of American Pathologists20	Hardy Diagnostics.....21
Cytomegalovirus (CMV) DNA by PCR, quantitative	Digoxin, assays and therapeutic drug monitoring (DIG)	Drugs of abuse detection	Enterovirus detection by RT-PCR
DiaSorin Molecular.....20	Randox Laboratories-US, Ltd.....24	College of American Pathologists20	DiaSorin Molecular.....20
Seegene Technologies.....24	Digoxin, serum	Randox Laboratories-US, Ltd.....24	Luminex.....22
Cytomegalovirus (CMV) DNA detection and quantitation	Carolina Liquid Chemistries Corp.....19	Drugs of abuse detection, controls	Seegene Technologies.....24
DiaSorin Molecular.....20	Randox Laboratories-US, Ltd.....24	Randox Laboratories-US, Ltd.....24	Enterovirus testing
Cytomegalovirus (CMV) DNA PCR, qualitative	Direct bilirubin (DBIL)	Drugs of abuse detection, immunoassays	DiaSorin Molecular.....20
Seegene Technologies.....24	Carolina Liquid Chemistries Corp.....19	Randox Laboratories-US, Ltd.....24	Luminex.....22
Cytomegalovirus (CMV) testing	Direct HDL cholesterol	Drugs of abuse detection, on-site kits	Seegene Technologies.....24
Meridian Bioscience Inc22	Randox Laboratories-US, Ltd.....24	Randox Laboratories-US, Ltd.....24	Enterovirus, real-time PCR
Cytopathology	Direct LDL cholesterol	Drugs of abuse detection, urine	DiaSorin Molecular.....20
College of American Pathologists20	Randox Laboratories-US, Ltd.....24	Randox Laboratories-US, Ltd.....24	Luminex.....22
	Disposable lab coats	Drugs of abuse detection, whole blood	Seegene Technologies.....24
	Hardy Diagnostics.....21	Randox Laboratories-US, Ltd.....24	Enzyme immunoassays
	Health Care Logistics21		IMMY21
	Disposable laboratory products		Randox Laboratories-US, Ltd.....24
	Hardy Diagnostics.....21		Enzyme immunoassays, automated
	Health Care Logistics21		Randox Laboratories-US, Ltd.....24
	DNA		Enzyme immunoassays, controls
	Luminex.....22		Randox Laboratories-US, Ltd.....24
	DNA ploidy analysis		Enzyme immunoassays, instrumentation
	College of American Pathologists20		Randox Laboratories-US, Ltd.....24
	DNA, amplification products		Enzymes, research use
	DiaSorin Molecular.....20		Randox Laboratories-US, Ltd.....24
	Seegene Technologies.....24		Epidermal growth factor receptor (EGFR) pharmDx™ (HER-1) oncoprotein, paraffin-embedded colorectal tissue
	DNA, extraction kits		Randox Laboratories-US, Ltd.....24
	Seegene Technologies.....24		Epidermal growth factor receptor (EGFR), paraffin block
	DNA, reagent kit		Randox Laboratories-US, Ltd.....24
	DiaSorin Molecular.....20		Epstein-Barr virus (EBV), DNA probe
	Seegene Technologies.....24		DiaSorin Molecular.....20
	DNA/RNA based, reagents		Error prevention
	DiaSorin Molecular.....20		Clinical Software Solutions.....19
	Meridian Bioscience Inc22		College of American Pathologists20
	Randox Laboratories-US, Ltd.....24		Escherichia coli
	Seegene Technologies.....24		Hardy Diagnostics.....21
	DNA/RNA based, systems		Luminex.....22
	Puritan Medical Products23		Seegene Technologies.....24
	Seegene Technologies.....24		Escherichia coli, 0157:H7
	Drug abuse profile (HHS/DOT), urine (five drugs)		Luminex.....22
	Randox Laboratories-US, Ltd.....24		Seegene Technologies.....24
	Drug abuse profile (HHS/DOT), urine (five drugs) (Includes GC/MS confirmation)		Escherichia coli, shiga toxin
	Randox Laboratories-US, Ltd.....24		Luminex.....22
	Drug abuse profile (routine), urine (eight drugs)		Esoteric testing
	Randox Laboratories-US, Ltd.....24		College of American Pathologists20
	Drug abuse profile (routine), urine (eight drugs) (Includes GC/MS confirmation)		Estradiol
	Randox Laboratories-US, Ltd.....24		Beckman Coulter, Inc.19
	Drug abuse profile, urine (eight drugs and alcohol)		Randox Laboratories-US, Ltd.....24
	Randox Laboratories-US, Ltd.....24		Estradiol assays
	Drug profile, blood (five drugs)		Beckman Coulter, Inc.19
	Randox Laboratories-US, Ltd.....24		Randox Laboratories-US, Ltd.....24
	Drug profile, blood (seven drugs)		Ethanol assays (See Alcohol assays)
	Randox Laboratories-US, Ltd.....24		Sekisui Diagnostics.....24
	Drug profile, qualitative, urine or gastric contents (with cannabinoids)		
	Randox Laboratories-US, Ltd.....24		

Page	Page	Page	Page
F			
Factor II (prothrombin) mutation analysis Seegene Technologies.....24	Fungus culture with stain Hardy Diagnostics.....21	Growth factor assays Randox Laboratories-US, Ltd.....24	Herpes simplex virus (HSV) by DNA probe DiaSorin Molecular.....20
Factor V Leiden mutation analysis DiaSorin Molecular.....20	Fungus stain Hardy Diagnostics.....21	Gynecologic pap test, liquid-based preparation and chlamydia trachomatis by nucleic acid amplification Seegene Technologies.....24	Herpes simplex virus (HSV) by DNA probe and quant PCR Seegene Technologies.....24
Factor V, Leiden mutation DNA test DiaSorin Molecular.....20 Seegene Technologies.....24	Furniture Health Care Logistics.....21 Proline Laboratory Furniture.....23	Gynecologic pap testing College of American Pathologists.....20	Herpes simplex virus (HSV) by PCR probe and quant PCR DiaSorin Molecular.....20 Seegene Technologies.....24
Fatty acids, free (nonester) Randox Laboratories-US, Ltd.....24	G	H	Herpes simplex virus (HSV) types I/II, DNA, PCR DiaSorin Molecular.....20 Luminex.....22 Seegene Technologies.....24
Fecal occult blood (FOB) (See Stool, occult blood) Hardy Diagnostics.....21	GAD65 (See Glutamic acid decarboxylase antibody kits) KRONUS, Inc.....21	H. pylori, serum or plasma Sekisui Diagnostics.....24	Herpes simplex virus testing Luminex.....22 Meridian Bioscience Inc.....22 Seegene Technologies.....24
Fecal occult blood (iFOB) Carolina Liquid Chemistries Corp.....19 Sekisui Diagnostics.....24	Gamma-glutamyl transpeptidase (GGT) Sekisui Diagnostics.....24	H. pylori, whole blood Sekisui Diagnostics.....24	Herpes virus testing Luminex.....22
Ferritin assay (FERR) Randox Laboratories-US, Ltd.....24	Genetic testing College of American Pathologists.....20	Haemophilus influenzae testing Seegene Technologies.....24	HerpeSelect, type-specific HSV antibody DiaSorin Molecular.....20
Fibrinogen Randox Laboratories-US, Ltd.....24	Genotyping Software SoftGenetics.....24	Hand scrub Health Care Logistics.....21	HHV6 DiaSorin Molecular.....20
Fibrinogen, assays Randox Laboratories-US, Ltd.....24	Gentamicin assay (GENT) Carolina Liquid Chemistries Corp.....19 Randox Laboratories-US, Ltd.....24	Hantavirus testing DiaSorin Molecular.....20	High density lipoprotein (HDL) cholesterol Randox Laboratories-US, Ltd.....24
Fine-needle aspiration cytology College of American Pathologists.....20	GGT Carolina Liquid Chemistries Corp.....19	Haptoglobin assays (HAPT) Carolina Liquid Chemistries Corp.....19	High sensitivity C-reactive protein (hs-CRP) Randox Laboratories-US, Ltd.....24
Flow cytometry College of American Pathologists.....20	Giardia lamblia detection Luminex.....22 Seegene Technologies.....24	Hazardous material storage Health Care Logistics.....21	High-density lipoprotein Randox Laboratories-US, Ltd.....24 Sekisui Diagnostics.....24
Folate (Folic acid), serum Randox Laboratories-US, Ltd.....24	Glassware Health Care Logistics.....21	HbA1c Beckman Coulter, Inc.....19	High-density lipoprotein, control Randox Laboratories-US, Ltd.....24
Folate assays Randox Laboratories-US, Ltd.....24	Glassware scrubbers Health Care Logistics.....21	HbA1c2 Beckman Coulter, Inc.....19	High-density lipoprotein, direct Randox Laboratories-US, Ltd.....24
Follicle-stimulating hormone (FSH) and luteinizing hormone (LH) Randox Laboratories-US, Ltd.....24	Gloves Health Care Logistics.....21	hCG (Human chorionic gonadotropin) Meridian Bioscience Inc.....22 Sekisui Diagnostics.....24	High-density lipoprotein, serum Randox Laboratories-US, Ltd.....24
Follicle-stimulating hormone assays Randox Laboratories-US, Ltd.....24	Glucose (GLU) Carolina Liquid Chemistries Corp.....19 Randox Laboratories-US, Ltd.....24 Sekisui Diagnostics.....24	hCG Ultra Combo Sekisui Diagnostics.....24	High-density lipoprotein, whole blood Randox Laboratories-US, Ltd.....24
Forensic toxicology College of American Pathologists.....20 Randox Laboratories-US, Ltd.....24	Glucose monitoring Randox Laboratories-US, Ltd.....24	HDL Carolina Liquid Chemistries Corp.....19	HIPAA compliance Accreditation Commission for Health Care (ACHC).....18 Clinical Software Solutions.....19
Fragile X Analysis Software SoftGenetics.....24	Glucose, assays Carolina Liquid Chemistries Corp.....19 Randox Laboratories-US, Ltd.....24	Helicobacter pylori detection Meridian Bioscience Inc.....22	Histocompatibility testing College of American Pathologists.....20
Free fatty acids Randox Laboratories-US, Ltd.....24	Glucose, body fluid Carolina Liquid Chemistries Corp.....19	Helicobacter pylori stool antigen Meridian Bioscience Inc.....22	Histology College of American Pathologists.....20
Free T3 assays Randox Laboratories-US, Ltd.....24	Glucose, standard Randox Laboratories-US, Ltd.....24	Helicobacter pylori urea breath test Meridian Bioscience Inc.....22	Histology Special Stains, Microorganisms Hardy Diagnostics.....21
Free T4 assays Randox Laboratories-US, Ltd.....24	Glutamic acid decarboxylase (GAD) autoantibody KRONUS, Inc.....21	Hematology products and services Beckman Coulter, Inc.....19 College of American Pathologists.....20 Electronic Imaging Materials.....20 Hardy Diagnostics.....21	Histopathology College of American Pathologists.....20
Free T4/TSH Randox Laboratories-US, Ltd.....24	Glutamic acid decarboxylase antibody kits KRONUS, Inc.....21	Hematology products and services, analyzers Beckman Coulter, Inc.....19 Carolina Liquid Chemistries Corp.....19	Histoplasma antibodies MiraVista Diagnostics.....23
Freezers Health Care Logistics.....21	Gram slide stainers Hardy Diagnostics.....21	Hematology products and services, controls American Proficiency Institute.....19 Beckman Coulter, Inc.....19	Histoplasma Antibody by Immunodiffusion IMMY.....21 MiraVista Diagnostics.....23
Fructosamine Carolina Liquid Chemistries Corp.....19	Gram's stain Hardy Diagnostics.....21	Hematopathology College of American Pathologists.....20	Histoplasma antigen detection IMMY.....21 MiraVista Diagnostics.....23
FT3 Randox Laboratories-US, Ltd.....24	Group A streptococcus Hardy Diagnostics.....21 Luminex.....22 Meridian Bioscience Inc.....22	Hemoglobin (Hgb) A1c Beckman Coulter, Inc.....19 Carolina Liquid Chemistries Corp.....19 Randox Laboratories-US, Ltd.....24	Histoplasma antigen test IMMY.....21 MiraVista Diagnostics.....23
FT4 Randox Laboratories-US, Ltd.....24	Group A streptococcus direct test, DNA probe Luminex.....22	Hemoglobin A1c Beckman Coulter, Inc.....19 Randox Laboratories-US, Ltd.....24	Histoplasma capsulatum IgG antibody MiraVista Diagnostics.....23
Fungal antibodies IMMY.....21	Group B strep PCR Luminex.....22 Meridian Bioscience Inc.....22 Seegene Technologies.....24	Hemoglobin HbA,c Randox Laboratories-US, Ltd.....24	
Fungal antibodies, qualitative IMMY.....21	Group B streptococcus colonization detection culture Hardy Diagnostics.....21	Hemoglobin, glycated assay Randox Laboratories-US, Ltd.....24	
Fungal antibodies, quantitative IMMY.....21			
Fungal serology College of American Pathologists.....20 IMMY.....21			
Fungus (mycology) culture Hardy Diagnostics.....21			

Page	Page	Page	Page
Homocyst(e)ine, plasma or serum Randox Laboratories-US, Ltd.24	Immunoassays, rapid Beckman Coulter, Inc.19 IMMY21	Influenza A and B, Real-time RT-PCR With Reflex to Influenza A, H1N1 (Novel) Luminex22 Seegene Technologies24	Randox Laboratories-US, Ltd.24
Homocysteine assay Randox Laboratories-US, Ltd.24	Immunoassays, reagents Beckman Coulter, Inc.19 IMMY21 Meridian Bioscience Inc22 Randox Laboratories-US, Ltd.24	Influenza A&B Control Kit, Rapid Sekisui Diagnostics24	Internet-based ordering and reporting Clinical Software Solutions19 College of American Pathologists20
hsCRP Carolina Liquid Chemistries Corp.19	Immunoblot kits DiaSorin Molecular20	Influenza A&B Kit, Rapid Meridian Bioscience Inc22	Iron Carolina Liquid Chemistries Corp.19 Sekisui Diagnostics24
hsTropoin Assay Beckman Coulter, Inc.19	Immunochemistry DiaSorin Molecular20 Randox Laboratories-US, Ltd.24	Influenza A, H1N1 (Novel), Real-time RT-PCR Seegene Technologies24	Iron assays (IRON) Randox Laboratories-US, Ltd.24
Human Identity Software SoftGenetics24	Immunodiagnosics Randox Laboratories-US, Ltd.24	Influenza Control Kit Microbiologics22 Sekisui Diagnostics24	Iron-binding capacity Carolina Liquid Chemistries Corp.19 Randox Laboratories-US, Ltd.24
Human Metapneumovirus Seegene Technologies24	Immunodiffusion kits IMMY21	Influenza Flu Test Kit Luminex22 Sekisui Diagnostics24	ISO 15189 Medical Laboratory Accreditation Services College of American Pathologists20
Human metapneumovirus (hMPV) RNA, PCR Seegene Technologies24	Immunoglobulin G (IGG) Beckman Coulter, Inc.19	Influenza Nasal Wash Kit Sekisui Diagnostics24	
Human papillomavirus (HPV) by PCR Seegene Technologies24	Immunoglobulin G, quantitative, cerebrospinal fluid (CSF) Beckman Coulter, Inc.19	Influenza syndrome panel Seegene Technologies24	Ketone bodies, serum Randox Laboratories-US, Ltd.24
Human papillomavirus (HPV) RNA assay Seegene Technologies24	Immunoglobulin G, quantitative, serum Beckman Coulter, Inc.19	Influenza testing DiaSorin Molecular20 Luminex22 Sekisui Diagnostics24	Kidney function, assessment Randox Laboratories-US, Ltd.24
I	Immunoglobulin G, specific Beckman Coulter, Inc.19	Influenza Virus Types A and B, RSV, Real-time PCR Luminex22	Kits, enzyme immunoassays IMMY21 Randox Laboratories-US, Ltd.24
IA-2 antibody kits KRONUS, Inc.21	Immunoglobulin M, quantitative, serum Beckman Coulter, Inc.19	Information management Advanced Data Systems Corporation18 Beckman Coulter, Inc.19 Clinical Software Solutions19 CompuGroup Medical US20 Gestalt20	L
IA-2 autoantibodies KRONUS, Inc.21	Immunoglobulins Beckman Coulter, Inc.19	Information systems Advanced Data Systems Corporation18 Beckman Coulter, Inc.19 Clinical Software Solutions19 CompuGroup Medical US20	Lab Benches Health Care Logistics21 Proline Laboratory Furniture23
Identity products and services Electronic Imaging Materials20	Immunoglobulins, assay Beckman Coulter, Inc.19	Insulin Carolina Liquid Chemistries Corp.19	Lab Carts Health Care Logistics21
iFOBT Controls Sekisui Diagnostics24	Immunohematology products Electronic Imaging Materials20	Insulin autoantibodies (IAA) KRONUS, Inc.21	Lab Coats Hardy Diagnostics21 Health Care Logistics21
iFOBT Patient Mailers Sekisui Diagnostics24	Immunohistochemistry College of American Pathologists20	Insulin, autoantibody kits KRONUS, Inc.21	Lab Furniture Health Care Logistics21
iFOBT Stool Collection Tubes Sekisui Diagnostics24	Immunology products and services American Proficiency Institute19 College of American Pathologists20 Electronic Imaging Materials20 Randox Laboratories-US, Ltd.24	Insulin-like growth factor I (IGF-I) Randox Laboratories-US, Ltd.24	Lab portals Clinical Software Solutions19 CompuGroup Medical US20
IgA Carolina Liquid Chemistries Corp.19 Randox Laboratories-US, Ltd.24	Immunology products and services, controls Randox Laboratories-US, Ltd.24	Insulin-like growth factor II (IGF-II) Randox Laboratories-US, Ltd.24	Labels CompuGroup Medical US20 Electronic Imaging Materials20 Health Care Logistics21
IgE Randox Laboratories-US, Ltd.24	Immunopathology products and services Electronic Imaging Materials20	Insulin-like Growth Factor-binding Protein 2 (IGFBP-2) Randox Laboratories-US, Ltd.24	Laboratory accreditation Accreditation Commission for Health Care (ACHC)18 College of American Pathologists20
IgG Beckman Coulter, Inc.19 Carolina Liquid Chemistries Corp.19 Hardy Diagnostics21 Randox Laboratories-US, Ltd.24	Incubators Health Care Logistics21	Interferon beta, antibody test Randox Laboratories-US, Ltd.24	Laboratory Accreditation Services Accreditation Commission for Health Care (ACHC)18 College of American Pathologists20
IgM Beckman Coulter, Inc.19 Carolina Liquid Chemistries Corp.19	Individualized Quality Control Plan (IQCP) Implementation Tools Clinical Software Solutions19 College of American Pathologists20	Interleukin beta, antibody test Randox Laboratories-US, Ltd.24	Laboratory automation (See also Automated analyzers) Advanced Data Systems Corporation18 Seegene Technologies24
Image Analysis Gestalt20	Infection control Hardy Diagnostics21	Interleukin-12 (IL-12) Randox Laboratories-US, Ltd.24	Laboratory automation/robotics Beckman Coulter, Inc.19 Seegene Technologies24
Image archival and storage Gestalt20	Infectious disease College of American Pathologists20 IMMY21 Randox Laboratories-US, Ltd.24	Interleukin-2 (IL-2) Randox Laboratories-US, Ltd.24	Laboratory carts Health Care Logistics21
Image Management Systems Gestalt20	Infectious disease, testing DiaSorin Molecular20 IMMY21	Interleukin-2 (IL-2) Randox Laboratories-US, Ltd.24	Laboratory information systems American Proficiency Institute19 Clinical Software Solutions19 CompuGroup Medical US20 Gestalt20 Orchard Software Corporation23
Immunoassays Beckman Coulter, Inc.19 IMMY21 Randox Laboratories-US, Ltd.24	Influenza A and B antibodies, quantitative Sekisui Diagnostics24	Interleukin-2 receptor (IL-2R) Randox Laboratories-US, Ltd.24	Laboratory process analysis Beckman Coulter, Inc.19 Clinical Software Solutions19
Immunoassays, analyzers Beckman Coulter, Inc.19 Carolina Liquid Chemistries Corp.19 Randox Laboratories-US, Ltd.24	Influenza A and B by direct immunoassay Sekisui Diagnostics24	Interleukin-6 Beckman Coulter, Inc.19 Randox Laboratories-US, Ltd.24	
Immunoassays, automated systems Beckman Coulter, Inc.19 Randox Laboratories-US, Ltd.24		Interleukin-6 (IL-6) Beckman Coulter, Inc.19	
Immunoassays, controls Randox Laboratories-US, Ltd.24			
Immunoassays, instrumentation Beckman Coulter, Inc.19 Randox Laboratories-US, Ltd.24			

Page	Page
College of American Pathologists20	
Gestalt.....20	
Laboratory safety	
College of American Pathologists20	
Health Care Logistics21	
Laboratory supplies	
Electronic Imaging Materials20	
Hardy Diagnostics.....21	
Health Care Logistics21	
Randox Laboratories-US, Ltd.....24	
Lactate dehydrogenase assays (IFCC) (LDH)	
Carolina Liquid Chemistries Corp.....19	
Sekisui Diagnostics.....24	
Latex agglutination tests	
Hardy Diagnostics.....21	
IMMY21	
LDH (Lactate dehydrogenase)	
Carolina Liquid Chemistries Corp.....19	
LDL	
Carolina Liquid Chemistries Corp.....19	
Lead, whole blood	
Meridian Bioscience Inc22	
Legionella detection	
Hardy Diagnostics.....21	
Luminex22	
Meridian Bioscience Inc22	
Seegene Technologies.....24	
Legionella pneumophila urine antigen	
Hardy Diagnostics.....21	
Leukocyte Adhesion Profile (CD11a, Cd11b, CD11c, CD18, CD45)	
Randox Laboratories-US, Ltd.....24	
Leukocytes	
Randox Laboratories-US, Ltd.....24	
Lidocaine assays	
Randox Laboratories-US, Ltd.....24	
Linearity materials	
Bio-Rad Laboratories19	
College of American Pathologists20	
Randox Laboratories-US, Ltd.....24	
Lipase	
Carolina Liquid Chemistries Corp.....19	
Sekisui Diagnostics.....24	
Lipase assay (LIP)	
Randox Laboratories-US, Ltd.....24	
Lipase, serum (LIP)	
Carolina Liquid Chemistries Corp.....19	
Randox Laboratories-US, Ltd.....24	
Lipid panel	
Randox Laboratories-US, Ltd.....24	
Lipid panel with LDL/HDL ratio	
Randox Laboratories-US, Ltd.....24	
Lipid panel with total cholesterol/HDL ratio	
Randox Laboratories-US, Ltd.....24	
Lipid testing	
Randox Laboratories-US, Ltd.....24	
Lipid testing, controls	
Randox Laboratories-US, Ltd.....24	
Lipids, total, serum	
Randox Laboratories-US, Ltd.....24	
Lipoprotein (a) (See LP(a))	
Carolina Liquid Chemistries Corp.....19	
Randox Laboratories-US, Ltd.....24	
Lipoproteins, assays	
Randox Laboratories-US, Ltd.....24	
Lithium (LI)	
Carolina Liquid Chemistries Corp.....19	
Lithium, serum	
Carolina Liquid Chemistries Corp.....19	

Page	Page
Randox Laboratories-US, Ltd.....24	
Low-density lipoprotein cholesterol (direct)	
Randox Laboratories-US, Ltd.....24	
Sekisui Diagnostics.....24	
Low-density lipoprotein testing	
Randox Laboratories-US, Ltd.....24	
Sekisui Diagnostics.....24	
LP(a)	
Carolina Liquid Chemistries Corp.....19	
Randox Laboratories-US, Ltd.....24	
M	
Magnesium	
Carolina Liquid Chemistries Corp.....19	
Sekisui Diagnostics.....24	
Magnesium, serum	
Carolina Liquid Chemistries Corp.....19	
Marijuana, urine assays (Cannabinoids)	
Randox Laboratories-US, Ltd.....24	
Maternal serum, control	
Randox Laboratories-US, Ltd.....24	
Media, chromogenic	
Hardy Diagnostics.....21	
Methadone assays	
Carolina Liquid Chemistries Corp.....19	
Randox Laboratories-US, Ltd.....24	
Methadone confirmation, urine	
Randox Laboratories-US, Ltd.....24	
Methadone screen, urine	
Carolina Liquid Chemistries Corp.....19	
Methamphetamine assays	
Randox Laboratories-US, Ltd.....24	
Methicillin-resistant staphylococcus aureus (MRSA) culture only Rare	
Hardy Diagnostics.....21	
Methicillin-resistant Staphylococcus aureus (MRSA), Nucleic Acid Amplification (NAA)	
Luminex22	
Seegene Technologies.....24	
Method development and validation	
Hardy Diagnostics.....21	
Microalbumin assays (MALB)	
Carolina Liquid Chemistries Corp.....19	
Randox Laboratories-US, Ltd.....24	
Sekisui Diagnostics.....24	
Microalbumin assays, controls	
Randox Laboratories-US, Ltd.....24	
Microbes, susceptibility testing (See Drug susceptibility testing)	
Hardy Diagnostics.....21	
Microbial strain typing	
Hardy Diagnostics.....21	
Microbiology products and services	
Beckman Coulter, Inc.19	
College of American Pathologists20	
Electronic Imaging Materials20	
Hardy Diagnostics.....21	
Microbiologics.....22	
Puritan Medical Products23	
Microbiology products and services, automated	
American Proficiency Institute.....19	
Beckman Coulter, Inc.19	
Electronic Imaging Materials20	
Hardy Diagnostics.....21	

Page	Page
Microorganisms	
Microbiologics.....22	
Microplate processing, automated	
Beckman Coulter, Inc.19	
Microplate processing, automated, calibrators	
Beckman Coulter, Inc.19	
Microsatellite Instability Software	
SoftGenetics.....24	
MLPA Software	
SoftGenetics.....24	
Molecular biology, products & services	
Electronic Imaging Materials20	
Microbiologics22	
Puritan Medical Products23	
Molecular biology, quality control	
Microbiologics22	
Molecular cytogenetics	
College of American Pathologists20	
Molecular diagnostics	
College of American Pathologists20	
DiaSorin Molecular.....20	
Hardy Diagnostics.....21	
Luminex22	
Randox Laboratories-US, Ltd.....24	
Molecular genetics	
College of American Pathologists20	
Luminex22	
Molecular pathology	
College of American Pathologists20	
Molecular pathology, controls	
Seegene Technologies.....24	
Mononucleosis	
Sekisui Diagnostics.....24	
Mononucleosis, testing	
Meridian Bioscience Inc22	
Morphine assay	
Randox Laboratories-US, Ltd.....24	
MRSA (See Staphylococcus indetification, methicillin resistant)	
Luminex22	
MTHFR & Hemochromatosis	
Seegene Technologies.....24	
MTHFR mutation analysis	
Seegene Technologies.....24	
Mycobacteriology	
College of American Pathologists20	
Mycobacterium tuberculosis (MTB) detection by PCR	
Seegene Technologies.....24	
Mycobacterium tuberculosis testing	
Seegene Technologies.....24	
Mycology products and services	
College of American Pathologists20	
Electronic Imaging Materials20	
Mycoplasma pneumoniae	
Luminex22	
Meridian Bioscience Inc22	
Seegene Technologies.....24	
Mycoplasma pneumoniae culture	
Hardy Diagnostics.....21	
Mycoplasma testing	
Luminex22	
Seegene Technologies.....24	
Mycoplasma/ureaplasma real-time DNA	
Seegene Technologies.....24	
Myoglobin	
Randox Laboratories-US, Ltd.....24	

Page	Page
Myoglobin, assays	
Randox Laboratories-US, Ltd.....24	
Myoglobin, controls	
Randox Laboratories-US, Ltd.....24	
N	
Needles	
Health Care Logistics21	
Neisseria gonorrhoeae by DNA probe	
Seegene Technologies.....24	
Neisseria gonorrhoeae testing	
Seegene Technologies.....24	
Nephelometry, kits	
Beckman Coulter, Inc.19	
Randox Laboratories-US, Ltd.....24	
Nephrology products and services	
Electronic Imaging Materials20	
Neurology products and services	
Electronic Imaging Materials20	
Neuropathology	
College of American Pathologists20	
Next Generation Sequencing Software	
Orchard Software Corporation.....23	
SoftGenetics.....24	
Nongynecologic pap testing	
College of American Pathologists20	
Norovirus detection, RT-PCR	
Seegene Technologies.....24	
Norovirus-RNA, PCR	
Seegene Technologies.....24	
Nucleic acid amplification products and detection systems (See also Polymerase chain reaction)	
College of American Pathologists20	
DiaSorin Molecular.....20	
Luminex22	
O	
Occupational/environmental toxicology	
Randox Laboratories-US, Ltd.....24	
OEM manufacturing	
Microbiologics.....22	
Puritan Medical Products23	
Randox Laboratories-US, Ltd.....24	
Oncology, products and services	
Gestalt.....20	
Randox Laboratories-US, Ltd.....24	
Opiate confirmation, urine	
Randox Laboratories-US, Ltd.....24	
Opiate screen, urine	
Carolina Liquid Chemistries Corp.....19	
Randox Laboratories-US, Ltd.....24	
Opiates, assays	
Carolina Liquid Chemistries Corp.....19	
Randox Laboratories-US, Ltd.....24	
Opiates, blood	
Randox Laboratories-US, Ltd.....24	
Organism identification, anaerobic bacteria	
Hardy Diagnostics.....21	
Organism identification, bacteria	
Hardy Diagnostics.....21	
Organism identification, mold	
Hardy Diagnostics.....21	

Page	Page	Page	Page
Organism identification, mycobacteria Seegene Technologies.....24	PCP (See Phencyclidine) Carolina Liquid Chemistries Corp.....19	Procalcitonin (PCT) Assay Beckman Coulter, Inc.....19	Randox Laboratories-US, Ltd.....24
Organism identification, yeast Hardy Diagnostics.....21	Pharmacogenomics (PGx) Luminex.....22	Proficiency testing American Proficiency Institute.....19	Quality assurance/quality control, molecular biology Bio-Rad Laboratories.....19
Oxycodone screen, urine Carolina Liquid Chemistries Corp.....19	Pharmacy management information Clinical Software Solutions.....19	College of American Pathologists.....20	Clinical Software Solutions.....19
Oxycodone/Oxymorphone Screen, Urine 100 Carolina Liquid Chemistries Corp.....19	Phencyclidine (PCP), qualitative, urine Carolina Liquid Chemistries Corp.....19	WSLH Proficiency Testing.....24	Microbiologies.....22
Oxycodone/Oxymorphone Screen, Urine 300 Carolina Liquid Chemistries Corp.....19	Phenobarbital, assays and therapeutic drug monitoring (PBARB) Carolina Liquid Chemistries Corp.....19	Prolactin assays Randox Laboratories-US, Ltd.....24	Quality assurance/quality control, products and services Beckman Coulter, Inc.....19
P	Phenytoin, assays and therapeutic drug monitoring (PHENY) Carolina Liquid Chemistries Corp.....19	Propoxyphene screen, urine Carolina Liquid Chemistries Corp.....19	Bio-Rad Laboratories.....19
Pain Management and Reports CompuGroup Medical US.....20	Phlebotomy products and supplies Electronic Imaging Materials.....20	Prostate cancer profile II (serial monitor) Beckman Coulter, Inc.....19	Clinical Software Solutions.....19
Pancreatic Amylase, Serum Carolina Liquid Chemistries Corp.....19	Phosphorus Carolina Liquid Chemistries Corp.....19	Prostate-specific antigen Beckman Coulter, Inc.....19	Microbiologies.....22
Panels Beckman Coulter, Inc.....19	Phosphorus, assays Randox Laboratories-US, Ltd.....24	Prostate-specific antigen (PSA), complexed Beckman Coulter, Inc.....19	Randox Laboratories-US, Ltd.....24
Panels, custom Beckman Coulter, Inc.....19	Phosphorus, urine Carolina Liquid Chemistries Corp.....19	Prostate-specific antigen (PSA), free/total ratio Beckman Coulter, Inc.....19	WSLH Proficiency Testing.....24
Panels, quality control Beckman Coulter, Inc.....19	Pipette tips Hardy Diagnostics.....21	Prostate-specific antigen (PSA), free/total ratio (serial monitor) Beckman Coulter, Inc.....19	Quality assurance/quality control, software Beckman Coulter, Inc.....19
Bio-Rad Laboratories.....19	Health Care Logistics.....21	Prostate-specific antigen (PSA), free/total ratio (without serial monitor) Beckman Coulter, Inc.....19	Bio-Rad Laboratories.....19
CompuGroup Medical US.....20	Luminex.....22	Prostate-specific antigen (PSA), free/total ratio (without serial monitor) Beckman Coulter, Inc.....19	Clinical Software Solutions.....19
Microbiologies.....22	Pipettes Hardy Diagnostics.....21	Prostate-specific antigen (PSA), free/total ratio (without serial monitor) Beckman Coulter, Inc.....19	College of American Pathologists.....20
Randox Laboratories-US, Ltd.....24	Health Care Logistics.....21	Prostate-specific antigen (PSA), free/total ratio (without serial monitor) Beckman Coulter, Inc.....19	CompuGroup Medical US.....20
Pap test products College of American Pathologists.....20	Luminex.....22	Prostate-specific antigen (PSA), free/total ratio (without serial monitor) Beckman Coulter, Inc.....19	Orchard Software Corporation.....23
Electronic Imaging Materials.....20	Pipettors Hardy Diagnostics.....21	Prostate-specific antigen (PSA), free/total ratio (without serial monitor) Beckman Coulter, Inc.....19	Randox Laboratories-US, Ltd.....24
Parainfluenza DiaSorin Molecular.....20	Health Care Logistics.....21	Prostate-specific antigen (PSA), free/total ratio (without serial monitor) Beckman Coulter, Inc.....19	Quality Management Tools Beckman Coulter, Inc.....19
Seegene Technologies.....24	Luminex.....22	Prostate-specific antigen (PSA), free/total ratio (without serial monitor) Beckman Coulter, Inc.....19	Clinical Software Solutions.....19
Parasite identification College of American Pathologists.....20	Point-of-care testing College of American Pathologists.....20	Prostate-specific antigen (PSA), free/total ratio (without serial monitor) Beckman Coulter, Inc.....19	College of American Pathologists.....20
Parasitology products and services College of American Pathologists.....20	Polymerase chain reaction College of American Pathologists.....20	Prostate-specific antigen (PSA), free/total ratio (without serial monitor) Beckman Coulter, Inc.....19	CompuGroup Medical US.....20
Electronic Imaging Materials.....20	Luminex.....22	Prostate-specific antigen (PSA), free/total ratio (without serial monitor) Beckman Coulter, Inc.....19	WSLH Proficiency Testing.....24
Parvovirus testing DiaSorin Molecular.....20	Potassium testing (K) (Electrolytes, analysis) Carolina Liquid Chemistries Corp.....19	Prostate-specific antigen (PSA), free/total ratio (without serial monitor) Beckman Coulter, Inc.....19	Quality System Management Beckman Coulter, Inc.....19
Paternity Analysis Software SoftGenetics.....24	Potassium testing (See also Electrolytes, analysis), standard Randox Laboratories-US, Ltd.....24	Prostate-specific antigen (PSA), free/total ratio (without serial monitor) Beckman Coulter, Inc.....19	Clinical Software Solutions.....19
Pathology Products and Services College of American Pathologists.....20	Potassium, serum Carolina Liquid Chemistries Corp.....19	Prostate-specific antigen (PSA), free/total ratio (without serial monitor) Beckman Coulter, Inc.....19	College of American Pathologists.....20
Electronic Imaging Materials.....20	Potassium, urine Carolina Liquid Chemistries Corp.....19	Prostate-specific antigen (PSA), free/total ratio (without serial monitor) Beckman Coulter, Inc.....19	
Gestalt.....20	Prealbumin (PALB) Carolina Liquid Chemistries Corp.....19	Prostate-specific antigen (PSA), free/total ratio (without serial monitor) Beckman Coulter, Inc.....19	R
Randox Laboratories-US, Ltd.....24	Prealbumin, nephelometry kit Randox Laboratories-US, Ltd.....24	Prostate-specific antigen (PSA), free/total ratio (without serial monitor) Beckman Coulter, Inc.....19	Radiology management information system Clinical Software Solutions.....19
Pathology products and services, anatomic College of American Pathologists.....20	Pregnancy test, urine Carolina Liquid Chemistries Corp.....19	Protein, controls Randox Laboratories-US, Ltd.....24	Reagent reservoirs Hardy Diagnostics.....21
Pathology products and services, clinical College of American Pathologists.....20	Randox Laboratories-US, Ltd.....24	Protein, standard Randox Laboratories-US, Ltd.....24	Reagents Beckman Coulter, Inc.....19
Pathology products and services, equipment Randox Laboratories-US, Ltd.....24	Pregnancy testing Meridian Bioscience Inc.....22	Protein, total Carolina Liquid Chemistries Corp.....19	Hardy Diagnostics.....21
Pathology products and services, molecular College of American Pathologists.....20	Randox Laboratories-US, Ltd.....24	Publications Accreditation Commission for Health Care (ACHC).....18	IMMY.....21
Pathology Workflow Gestalt.....20	Primidone, assays and therapeutic drug monitoring Carolina Liquid Chemistries Corp.....19	College of American Pathologists.....20	Luminex.....22
Patient management information Advanced Data Systems Corporation.....18	Procalcitonin (PCT) Assay Beckman Coulter, Inc.....19		Meridian Bioscience Inc.....22
Clinical Software Solutions.....19	Procalcitonin (PCT) Assay Beckman Coulter, Inc.....19		Randox Laboratories-US, Ltd.....24
College of American Pathologists.....20	Procalcitonin (PCT) Assay Beckman Coulter, Inc.....19		Reagents - coagulation assay Hardy Diagnostics.....21

	Page
Reagents, instrument-specific, Hitachi	
Randox Laboratories-US, Ltd.....	24
Reagents, lipid testing	
Beckman Coulter, Inc.....	19
Randox Laboratories-US, Ltd.....	24
Reagents, lipid testing system	
Randox Laboratories-US, Ltd.....	24
Reagents, stable liquid	
Randox Laboratories-US, Ltd.....	24
Reagents, therapeutic drug monitoring	
Randox Laboratories-US, Ltd.....	24
Reagents, virology	
Randox Laboratories-US, Ltd.....	24
Red blood cell (RBC) antigens	
Beckman Coulter, Inc.....	19
Red blood cell (RBC) count	
Beckman Coulter, Inc.....	19
Reference books	
College of American Pathologists.....	20
Reference laboratory testing	
Beckman Coulter, Inc.....	19
Reference materials	
Luminex.....	22
Refrigerators	
Health Care Logistics.....	21
Remote diagnostics	
Beckman Coulter, Inc.....	19
Renal function panel	
Randox Laboratories-US, Ltd.....	24
Reproductive hormones	
College of American Pathologists.....	20
Respiratory syncytial virus by PCR	
DiaSorin Molecular.....	20
Luminex.....	22
Seegene Technologies.....	24
Respiratory syncytial virus testing	
Luminex.....	22
Meridian Bioscience Inc.....	22
Seegene Technologies.....	24
Respiratory Virus Pan by RT-PCR	
Luminex.....	22
Respiratory Virus Panel (RVP), PCR	
Luminex.....	22
Respiratory virus tests	
Luminex.....	22
Rheumatoid factor assays (RF)	
Carolina Liquid Chemistries Corp.....	19
Randox Laboratories-US, Ltd.....	24
Rickettsia detection	
Seegene Technologies.....	24
Rickettsia testing	
DiaSorin Molecular.....	20
Seegene Technologies.....	24
Rotavirus detection	
Luminex.....	22
Meridian Bioscience Inc.....	22
Seegene Technologies.....	24
RSV	
Luminex.....	22
RSV (Respiratory syncytial virus)	
Luminex.....	22

S

Safety products	
Electronic Imaging Materials.....	20
Hardy Diagnostics.....	21
Health Care Logistics.....	21
Safety Products, Disposable	
Hardy Diagnostics.....	21

Health Care Logistics	21
Safety Products, Lab Coats	
Hardy Diagnostics.....	21
Health Care Logistics.....	21
Salicylate, serum	
Carolina Liquid Chemistries Corp.....	19
Salicylates, assays	
Carolina Liquid Chemistries Corp.....	19
Sekisui Diagnostics.....	24
Salmonella identification	
Hardy Diagnostics.....	21
Luminex.....	22
Seegene Technologies.....	24
Sample management systems	
Clinical Software Solutions.....	19
SARS Associated Coronavirus by RT-PCR	
College of American Pathologists.....	20
Luminex.....	22
Seegene Technologies.....	24
SARS coronavirus RNA, PCR	
College of American Pathologists.....	20
Luminex.....	22
Seegene Technologies.....	24
SARS-CoV-2 - Antibody Test	
Beckman Coulter, Inc.....	19
College of American Pathologists.....	20
Hardy Diagnostics.....	21
Luminex.....	22
Seminars	
College of American Pathologists.....	20
Sequence Analysis Software	
SoftGenetics.....	24
Serology products and services	
Electronic Imaging Materials.....	20
Randox Laboratories-US, Ltd.....	24
Serum proteins	
Randox Laboratories-US, Ltd.....	24
Serum proteins, calibrators	
Randox Laboratories-US, Ltd.....	24
Serum proteins, controls	
Randox Laboratories-US, Ltd.....	24
Serum, specialty controls	
Randox Laboratories-US, Ltd.....	24
Services	
American Proficiency Institute.....	19
Clinical Software Solutions.....	19
Sharps collection systems	
Health Care Logistics.....	21
Shigella identification	
Hardy Diagnostics.....	21
Luminex.....	22
Seegene Technologies.....	24
Slide stainers	
Hardy Diagnostics.....	21
Sodium testing (NA) (Electrolytes, analysis)	
Carolina Liquid Chemistries Corp.....	19
Sodium, serum	
Carolina Liquid Chemistries Corp.....	19
Sodium, urine	
Carolina Liquid Chemistries Corp.....	19
Software connectivity	
Clinical Software Solutions.....	19
CompuGroup Medical US.....	20
Gestalt.....	20
Sorbitol Dehydrogenase (SDH)	
Sekisui Diagnostics.....	24
Specimen management systems	
Clinical Software Solutions.....	19
College of American Pathologists.....	20
Specimen, collection	
Health Care Logistics.....	21

Puritan Medical Products	23
Specimen, transport	
Health Care Logistics.....	21
Puritan Medical Products.....	23
Specimens	
Puritan Medical Products.....	23
Specimens, transport systems	
Puritan Medical Products.....	23
Sporothrix antibody	
IMMY.....	21
Staffing	
Gestalt.....	20
Stains	
Hardy Diagnostics.....	21
Standards	
Randox Laboratories-US, Ltd.....	24
Staph mec A	
Luminex.....	22
Staphylococcus identification, methicillin resistant	
Hardy Diagnostics.....	21
Sterile production service	
Puritan Medical Products.....	23
Stool, occult blood	
Hardy Diagnostics.....	21
Storage equipment	
Health Care Logistics.....	21
Strep A antigen	
Sekisui Diagnostics.....	24
Strep A, Ultra	
Sekisui Diagnostics.....	24
Streptococcus - group A	
Hardy Diagnostics.....	21
Luminex.....	22
Streptococcus (group B) culture, genital	
Hardy Diagnostics.....	21
Streptococcus identification	
Hardy Diagnostics.....	21
Streptococcus pneumoniae antigen	
Hardy Diagnostics.....	21
Streptococcus pneumoniae testing	
Hardy Diagnostics.....	21
Luminex.....	22
Seegene Technologies.....	24
Surgical pathology	
College of American Pathologists.....	20
Susceptibility testing	
College of American Pathologists.....	20
Hardy Diagnostics.....	21
Syphilis, diagnosis	
Randox Laboratories-US, Ltd.....	24

T

T3	
Randox Laboratories-US, Ltd.....	24
T4	
Carolina Liquid Chemistries Corp.....	19
Randox Laboratories-US, Ltd.....	24
TDM (See Therapeutic drug monitoring)	
College of American Pathologists.....	20
Randox Laboratories-US, Ltd.....	24
Temperature Monitoring, Automated	
Health Care Logistics.....	21
Temperature Monitoring, Continuous	
Health Care Logistics.....	21

Temporary staffing	
Gestalt.....	20
Test kits, rapid	
Beckman Coulter, Inc.....	19
Hardy Diagnostics.....	21
IMMY.....	21
Test Tube Racks	
Health Care Logistics.....	21
IMMY.....	21
Testosterone, assays	
Carolina Liquid Chemistries Corp.....	19
Theophylline, assay and therapeutic drug monitoring (THEO)	
Carolina Liquid Chemistries Corp.....	19
Theophylline, serum	
Carolina Liquid Chemistries Corp.....	19
Therapeutic drug monitoring	
Carolina Liquid Chemistries Corp.....	19
College of American Pathologists.....	20
Therapeutic drug monitoring, controls	
Randox Laboratories-US, Ltd.....	24
Therapeutic drug monitoring, reagents	
Randox Laboratories-US, Ltd.....	24
Thyroglobulin, test kits	
Beckman Coulter, Inc.....	19
Randox Laboratories-US, Ltd.....	24
Thyroid function	
Randox Laboratories-US, Ltd.....	24
Thyroid hormone assays	
Carolina Liquid Chemistries Corp.....	19
Thyroid panel	
Randox Laboratories-US, Ltd.....	24
Thyroid panel with TSH	
Randox Laboratories-US, Ltd.....	24
Thyroid peroxidase, test kits	
Randox Laboratories-US, Ltd.....	24
Thyroid-stimulating hormone (TSH)	
Carolina Liquid Chemistries Corp.....	19
Randox Laboratories-US, Ltd.....	24
Thyroid-stimulating hormone receptor antibody tests	
KRONUS, Inc.....	21
Thyroxine (See also Thyroid hormone assays)	
Carolina Liquid Chemistries Corp.....	19
Thyroxine (T4)	
Carolina Liquid Chemistries Corp.....	19
TIBC (See Iron-binding capacity)	
Carolina Liquid Chemistries Corp.....	19
Tobramycin, assay and therapeutic drug monitoring (TOB)	
Carolina Liquid Chemistries Corp.....	19
Total IgE	
Randox Laboratories-US, Ltd.....	24
Toxicology products and services	
Clinical Software Solutions.....	19
College of American Pathologists.....	20
CompuGroup Medical US.....	20
Gestalt.....	20
Randox Laboratories-US, Ltd.....	24
Training	
Accreditation Commission for Health Care (ACHC).....	18
Clinical Software Solutions.....	19
WSLH Proficiency Testing.....	24
Transferrin (TRF)	
Carolina Liquid Chemistries Corp.....	19
Randox Laboratories-US, Ltd.....	24
Transferrin assays (TRF)	
Randox Laboratories-US, Ltd.....	24

	Page
Transfusion medicine	
College of American Pathologists	20
Treponema pallidum	
Seegene Technologies.....	24
Trichomonas vaginalis by Nucleic Acid Amplification (NAA)	
Seegene Technologies.....	24
Trichomonas vaginalis, rapid testing	
Sekisui Diagnostics.....	24
Tricyclic antidepressants, serum	
Carolina Liquid Chemistries Corp.....	19
Triglyceride (TRIG)	
Randox Laboratories-US, Ltd.....	24
Triglycerides	
Carolina Liquid Chemistries Corp.....	19
Randox Laboratories-US, Ltd.....	24
Sekisui Diagnostics.....	24
Triiodothyronine (T3)	
Randox Laboratories-US, Ltd.....	24
Triiodothyronine (T3), free, serum	
Randox Laboratories-US, Ltd.....	24
Trisomy/Aneuploidy Software	
SoftGenetics.....	24
Troponin I	
Beckman Coulter, Inc.	19
Randox Laboratories-US, Ltd.....	24
Troponin testing	
Randox Laboratories-US, Ltd.....	24
TSH (See Thyroid-stimulating hormone)	
Randox Laboratories-US, Ltd.....	24
Tubes	
Health Care Logistics	21
IMMY	21
Tumor markers	
College of American Pathologists	20
Randox Laboratories-US, Ltd.....	24
Tumor necrosis factor alfa	
Randox Laboratories-US, Ltd.....	24
T-uptake test (TUP)	
Carolina Liquid Chemistries Corp.....	19

U

Unsaturated Iron Binding Capacity (UIBC)	
Sekisui Diagnostics.....	24
Urea	
Randox Laboratories-US, Ltd.....	24
Sekisui Diagnostics.....	24
Ureaplasma testing	
Seegene Technologies.....	24
Ureaplasma/mycoplasma hominis culture	
Seegene Technologies.....	24
Uric acid (UA)	
Carolina Liquid Chemistries Corp.....	19
Randox Laboratories-US, Ltd.....	24
Sekisui Diagnostics.....	24
Uric acid, assays	
Beckman Coulter, Inc.	19
Uric acid, body fluid	
Beckman Coulter, Inc.	19
Uric acid, urine	
Beckman Coulter, Inc.	19
Urinalysis	
Beckman Coulter, Inc.	19
College of American Pathologists	20
Hardy Diagnostics.....	21
Randox Laboratories-US, Ltd.....	24

	Page
Urinalysis test reports	
Beckman Coulter, Inc.	19
Urinalysis, attended analyzer	
Beckman Coulter, Inc.	19
Urinalysis, automated	
Beckman Coulter, Inc.	19
Urinalysis, complete (with microscopic examination)	
Beckman Coulter, Inc.	19
Urinalysis, controls	
Beckman Coulter, Inc.	19
Randox Laboratories-US, Ltd.	24
Urinalysis, routine (with microscopic examination on positives)	
Beckman Coulter, Inc.	19
Urinalysis, specimen transport	
Beckman Coulter, Inc.	19
Urinary tract infections	
Beckman Coulter, Inc.	19
Urine	
Beckman Coulter, Inc.	19
Randox Laboratories-US, Ltd.	24
Urine chemistries, amylase	
Beckman Coulter, Inc.	19
Randox Laboratories-US, Ltd.	24
Urine chemistries, BUN	
Beckman Coulter, Inc.	19
Randox Laboratories-US, Ltd.	24
Urine chemistries, calcium	
Beckman Coulter, Inc.	19
Randox Laboratories-US, Ltd.	24
Urine chemistries, chloride	
Beckman Coulter, Inc.	19
Randox Laboratories-US, Ltd.	24
Urine chemistries, creatinine	
Beckman Coulter, Inc.	19
Randox Laboratories-US, Ltd.	24
Urine chemistries, glucose	
Beckman Coulter, Inc.	19
Randox Laboratories-US, Ltd.	24
Urine chemistries, magnesium	
Beckman Coulter, Inc.	19
Randox Laboratories-US, Ltd.	24
Urine chemistries, micro total protein	
Beckman Coulter, Inc.	19
Urine chemistries, pancreatic amylase	
Beckman Coulter, Inc.	19
Urine chemistries, phosphorus	
Beckman Coulter, Inc.	19
Randox Laboratories-US, Ltd.	24
Urine chemistries, potassium	
Beckman Coulter, Inc.	19
Randox Laboratories-US, Ltd.	24
Urine chemistries, sodium	
Beckman Coulter, Inc.	19
Randox Laboratories-US, Ltd.	24
Urine chemistries, urea	
Beckman Coulter, Inc.	19
Randox Laboratories-US, Ltd.	24
Urine chemistries, uric acid	
Beckman Coulter, Inc.	19
Randox Laboratories-US, Ltd.	24
Urine culture, comprehensive	
Beckman Coulter, Inc.	19
Urine culture, routine	
Beckman Coulter, Inc.	19
Urine cytology	
Beckman Coulter, Inc.	19

	Page
Urine protein, electrophoresis	
Beckman Coulter, Inc.	19
Urine, analyzers	
Beckman Coulter, Inc.	19
Randox Laboratories-US, Ltd.	24
Urine, cell analysis	
Beckman Coulter, Inc.	19
Urine, cell counting	
Beckman Coulter, Inc.	19
Urine, chemistry, controls	
Beckman Coulter, Inc.	19
Randox Laboratories-US, Ltd.	24
Urine, ELISA kits	
Randox Laboratories-US, Ltd.	24
Urine, human, control	
Randox Laboratories-US, Ltd.	24
Urine, nitrites	
Beckman Coulter, Inc.	19
Urine, oral fluid assays	
Beckman Coulter, Inc.	19
Urine, oxidants	
Beckman Coulter, Inc.	19
Urine, pH	
Beckman Coulter, Inc.	19
Urine, sediment analysis	
Beckman Coulter, Inc.	19
Urine, specialty controls	
Randox Laboratories-US, Ltd.	24
Urine, specific gravity	
Beckman Coulter, Inc.	19
Urine, specimen transport	
Health Care Logistics	21
Urine, test strips	
Beckman Coulter, Inc.	19

V

Vaccine Monitoring Solutions	
Gestalt.....	20
Valproic acid, assay and therapeutic drug monitoring (VPA)	
Carolina Liquid Chemistries Corp.....	19
VanA and VanB real-time PCR	
Seegene Technologies.....	24
Vancomycin assay (VANC)	
Carolina Liquid Chemistries Corp.....	19
Varicella-zoster virus testing	
DiaSorin Molecular.....	20
Luminex.....	22
Seegene Technologies.....	24
Viral load testing	
College of American Pathologists.....	20
Viral respiratory panel by PCR	
Seegene Technologies.....	24
Virology products and services	
College of American Pathologists.....	20
DiaSorin Molecular.....	20
Electronic Imaging Materials.....	20
Hardy Diagnostics.....	21
Microbiologies.....	22
Randox Laboratories-US, Ltd.....	24
Virus	
Seegene Technologies.....	24
Vitamin B1	
Beckman Coulter, Inc.....	19
Vitamin B1 (Thiamine)	
Beckman Coulter, Inc.....	19
Vitamin B1, plasma	
Beckman Coulter, Inc.....	19
Vitamin B1, whole blood	
Beckman Coulter, Inc.....	19

	Page
Vitamin B12	
Beckman Coulter, Inc.	19
Vitamin B12 : Holotranscobalamin Assess®	
Beckman Coulter, Inc.	19
Vitamin B12 and folate	
Beckman Coulter, Inc.	19
Vitamin B12 unsaturated binding capacity	
Beckman Coulter, Inc.	19
Vitamin B2	
Beckman Coulter, Inc.	19
Vitamin B2 (Riboflavin)	
Beckman Coulter, Inc.	19
Vitamin B6, plasma	
Beckman Coulter, Inc.	19
Vitamin D	
Beckman Coulter, Inc.	19
Vitamin D assay	
Beckman Coulter, Inc.	19
Vitamin D, 1,25 dihydroxy	
Beckman Coulter, Inc.	19
Vitamin D, 25-hydroxy	
Beckman Coulter, Inc.	19
Carolina Liquid Chemicals Corp.	19
Voltage-gated calcium channel (VGCC) antibody kits	
KRONUS, Inc.	21
VZV (See Varicella-zoster virus)	
Luminex	22

W

Waste Disposal	
Health Care Logistics	21
Water	
Health Care Logistics	21
WBC	
Beckman Coulter, Inc.	19
Web Based Courses	
Beckman Coulter, Inc.	19
College of American Pathologists	20
West Nile IgG & IgM, serum	
DiaSorin Molecular	20
West Nile virus	
DiaSorin Molecular	20
Western blot technology	
Randox Laboratories-US, Ltd.	24
Whole Slide Image Viewer	
Gestalt	20
Workshops	
Accreditation Commission for Health Care (ACHC)	18
Beckman Coulter, Inc.	19
College of American Pathologists	20
Luminex	22

Z

Zinc, serum	
Randox Laboratories-US, Ltd.....	24
Zinc, urine	
Randox Laboratories-US, Ltd.....	24

Individual EUAs for Serology and other Adaptive Immune Response Tests for SARS-CoV-2

Company Name	Last Updated (and Original Date EUA Issued)	Name of Test	Type of Test	Authorized Settings
Abbott Laboratories Inc.	5/24/2022, 4/26/2020	SARS-CoV-2 IgG assay	IgG, CMIA	H, M
Abbott Laboratories Inc.	4/7/2022, 3/1/2021	AdviseDx SARS-CoV-2 IgG II	IgG, CMIA, Semi-quantitative	H, M
Abbott Laboratories Inc.	10/1/2021	AdviseDx SARS-CoV-2 IgM 10/09/2020	IgM, CMIA	H, M
Access Bio, Inc.	8/12/2021, 7/24/2020	CareStart COVID-19 IgM/IgG	IgM and IgG, Lateral Flow, Fingerstick Whole Blood	H, M, W
Adaptive Biotechnologies Corporation	9/2/2021, 3/5/2021	T-Detect COVID Test	T-cell receptor beta (TCR β), Sequencing	H
Assure Tech. (Hangzhou Co., Ltd)	1/31/2022, 7/6/2020	Assure COVID-19 IgG/IgM Rapid Test Device	IgM and IgG, Lateral Flow, Fingerstick Whole Blood	H, M, W
Babson Diagnostics, Inc.	12/15/2021, 6/23/2020	Babson Diagnostics aC19G1	IgG, CLIA	H
Beckman Coulter, Inc.	12/17/2021, 6/26/2020	Access SARS-CoV-2 IgG	IgG, CLIA	H, M
Beckman Coulter, Inc.	12/13/2021, 10/8/2020	Access SARS-CoV-2 IgM	IgM, CLIA	H, M
Beckman Coulter, Inc.	8/18/2021, 3/3/2021	Access SARS-CoV-2 IgG II	IgG, CLIA, Semi-quantitative	H, M
Biohit Healthcare (Hefei) Co. Ltd.	12/16/2021, 6/18/2020	Biohit SARS-CoV-2 IgM/IgG Antibody Test Kit	IgM and IgG Lateral Flow	H, M
Bio-Rad Laboratories, Inc.	2/25/2022, 4/29/2020	Platelia SARS-CoV-2 Total Ab assay	Total Antibody, ELISA	H
Diabetomics, Inc.	8/10/2021, 6/4/2021	CovAb SARS-CoV-2 Ab Test	Total Antibody, Lateral Flow, Oral Fluid	H, M, W
DiaSorin Inc.	10/7/2021, 4/24/2020	LIAISON SARS-CoV-2 S1/S2 IgG	IgG, CLIA	H, M
Diazyme Laboratories, Inc.	1/18/2022, 7/8/2020	Diazyme DZ-Lite SARS-CoV-2 IgG CLIA Kit	IgG, CLIA, Semi-quantitative	H, M
Diazyme Laboratories, Inc.	1/18/2022, 8/17/2020	Diazyme DZ-Lite SARS-CoV-2 IgM CLIA Kit	IgM, CLIA	H, M
Emory Medical Laboratories	4/21/2022, 6/15/2020	SARS-CoV-2 RBD IgG test	IgG, ELISA	H
EUROIMMUN US Inc.	2/9/2022, 5/4/2020	Anti-SARS-CoV-2 ELISA (IgG)	IgG, ELISA	H
EUROIMMUN US, Inc.	3/1/2022, 10/4/2021	EUROIMMUN Anti-SARS-CoV-2 S1 Curve ELISA (IgG)	IgG, ELISA, Semi-quantitative	H, M
Hangzhou Laihe Biotech Co., Ltd.	2/14/2022, 6/19/2020	LYHER Novel Coronavirus (2019-nCoV) IgM/IgG Antibody Combo Test Kit (Colloidal Gold)	IgM and IgG Lateral Flow	H, M
Healgen Scientific LLC	1/26/2022, 5/29/2020	COVID-19 IgG/IgM Rapid Test Cassette (Whole Blood/Serum/Plasma)	IgM and IgG Lateral Flow	H, M
Immunodiagnostic Systems Ltd.	12/9/2021, 2/10/2021	IDS SARS-CoV-2 IgG	IgG, CLIA	H, M
InBios International, Inc	8/24/2021	SCoV-2 Detect IgG Rapid Test	IgG, Lateral Flow, Fingerstick Whole Blood	H, M, W
InBios International, Inc.	10/22/2021	SCoV-2 Detect Neutralizing Ab ELISA	Total Neutralizing Antibodies, ELISA	H
Jiangsu Well Biotech Co., Ltd.	3/1/2022, 9/30/2020	Orawell IgM/IgG Rapid Test	IgM and IgG, Lateral Flow	H, M
LG Chem, Ltd.	5/31/2022	AdvanSure SARS-CoV-2 IgG(RBD) ELISA	IgG, ELISA	H
LG Chem, Ltd.	5/19/2022	AdvanSure SARS-CoV-2 IgG(S1) ELISA	IgG, ELISA	H
Luminex Corporation	3/9/2022, 7/16/2020	xMAP SARS-CoV-2 Multi-Antigen IgG Assay	IgG, FMIA	H
LumiraDx UK Ltd.	8/2/2021	LumiraDx SARS-CoV-2 Ab Test	Total Antibody, Fluorescence Immunoassay, Fingerstick Whole Blood	H, M, W
Megna Health, Inc.	3/9/2022, 7/17/2020	Rapid COVID-19 IgM/IgG Combo Test Kit	IgM and IgG, Lateral Flow, Fingerstick Whole Blood	H, M, W
NanoEntek America, Inc.	12/13/2021, 9/29/2020	FREND COVID-19 total Ab	Total Antibody, FIA	H, M
Nirmidas Biotech, Inc.	11/2/2021, 12/31/2020	MidaSpot COVID-19 Antibody Combo Detection Kit	IgM and IgG, lateral flow, Fingerstick Whole Blood	H, M, W
NOWDiagnostics, Inc.	2/1/2022, 5/4/2021	ADEXUSdx COVID-19 Test	Total Antibody, Lateral Flow, Fingerstick Whole Blood	H, M, W
Ortho-Clinical Diagnostics, Inc.	9/8/2021, 7/9/2021	VITROS Immunodiagnostic Products Anti-SARS-CoV-2 IgG Quantitative Reagent Pack used in combination with the VITROS Immunodiagnostic Products Anti-SARS-CoV-2 IgG Quantitative Calibrator	IgG, CLIA, Quantitative	H, M
Ortho-Clinical Diagnostics, Inc.	7/22/2021	VITROS Immunodiagnostic Products Anti-SARS-CoV-2 Total N Reagent Pack used in combination with the VITROS Immunodiagnostic Products Anti-SARS-CoV-2 Total N Antibody Calibrators	Total Antibody, CLIA	H, M
QIAGEN, GmbH	2/28/2022, 5/11/2021	QIArearch Anti-SARS-CoV-2 Total Test	Total Antibody, Digital Lateral Flow	H, M

connect. engage.
discover.



The Lab Directors

SEPTEMBER 12-14, 2022 | AUSTIN, TEXAS
www.labdirectorssummit.com

S u m m i t



In partnership with



Save the date for 2023!
SEPTEMBER 20-22, 2023



Company Name	Last Updated (and Original Date EUA Issued)	Name of Test	Type of Test	Authorized Settings
Roche Diagnostics, Inc.	2/9/2022, 11/25/2020	Elecsys Anti-SARS-CoV-2 S	Total Antibody, ECLIA, Semi-quantitative	H, M
Shenzhen New Industries Biomedical Engineering Co., Ltd.	12/16/2021, 9/14/2020	MAGLUMI 2019-nCoV IgM/IgG	IgM and IgG, CLIA	H, M
Siemens Healthcare Diagnostics Inc.	3/4/2022, 1/8/2021	Dimension EXL SARS-CoV-2 IgG (CV2G)	IgG, CLIA, Semi-quantitative	H, M
Siemens Healthcare Diagnostics Inc.	3/4/2022, 1/8/2021	Dimension Vista SARS-CoV-2 IgG (COV2G)	IgG, CLIA, Semi-quantitative	H, M
Siemens Healthcare Diagnostics Inc.	1/19/2022, 3/23/2021	Atellica IM SARS-CoV-2 IgG (sCOVG)	IgG, CLIA, Semi-quantitative	H, M
Siemens Healthcare Diagnostics Inc.	1/18/2022, 6/8/2020	Dimension Vista SARS-CoV-2 Total antibody assay (COV2T)	Total Antibody, CLIA	H, M
Siemens Healthcare Diagnostics Inc.	1/18/2022, 6/8/2020	Dimension EXL SARS-CoV-2 Total antibody assay (CV2T)	Total Antibody, CLIA	H, M
Siemens Healthcare Diagnostics Inc.	9/14/2021, 5/29/2020	ADVIA Centaur SARS-CoV-2 Total (COV2T)	Total Antibody, CLIA, Semi-quantitative	H, M
Siemens Healthcare Diagnostics Inc.	8/10/2021, 6/17/2021	ADVIA Centaur SARS-CoV-2 IgG (sCOVG)	IgG, CLIA, Semi-quantitative	H, M
Sugentech, Inc.	11/19/2021, 9/3/2020	SGTi-flex COVID-19 IgG	IgG, Lateral Flow, Fingerstick Whole Blood	H, M, W
ZEUS Scientific, Inc.	12/13/2021, 5/11/2021	ZEUS ELISA SARS-CoV-2 Total Antibody Test System	Total Antibody, ELISA	H, M
ZEUS Scientific, Inc.	12/8/2021, 10/6/2020	ZEUS ELISA SARS-CoV-2 IgG Test System	IgG, ELISA	H, M

Authorized settings include the following:

H - Laboratories certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA), 42 U.S.C. §263a, that meet requirements to perform high complexity tests.

M - Laboratories certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA), 42 U.S.C. §263a, that meet requirements to perform moderate complexity tests.

W - Patient care settings operating under a CLIA Certificate of Waiver.

2 Authorization Documents include the Healthcare Provider (HCP) and Patient Fact Sheets and either the Manufacture Instructions/Package Insert (abbreviated to IFU) or the EUA Summary.

3 Abbreviations: CLIA = chemiluminescence immunoassay; ELISA = enzyme-linked immunosorbent assay; ECLIA = electrochemiluminescence immunoassay; FMIA = fluorescent microsphere immunoassay; CMIA = chemiluminescent microparticle immunoassay; ELFA = enzyme-linked fluorescence assay

Source: FDA: In Vitro Diagnostics EUAs - Serology and Other Adaptive Immune Response Tests for SARS-CoV-2

<https://www.fda.gov/medical-devices/coronavirus-disease-2019-covid-19-emergency-use-authorizations-medical-devices/in-vitro-diagnostics-euas-serology-and-other-adaptive-immune-response-tests-sars-cov-2> (Accessed June 29, 2022).

INDEX OF ADVERTISERS

ADVERTISER	PAGE
ACHC achc.org	21
American Proficiency Institute api-pt.com	29
College of American Pathologists cap.org	C2
CompuGroup Medical US..... cgm.com/us	23
Diasorin Molecular molecular.diasorin.com	C3
Health Care Logistics gohcl.com	40
KRONUS kronus.com.....	C4
LumiraDx - Fast Lab Solutions lumiradx.com	13
Randox Laboratories store.randox.com.....	3
Sekisui Diagnostics sekisuidiagnostics.com/crp..	15

This index is provided as a service.

The publisher does not assume liability for errors or omissions.



LOOKING FOR THE PERFECT PPE ORGANIZER?


LOOK NO FURTHER!

Health Care Logistics[®] INC.

GOHCL.COM • 1.800.848.1633

©Health Care Logistics, Inc. 2022

Luminex and DiaSorin Molecular Provide Molecular Diagnostic Solutions to Optimize Your Laboratory Services



Respiratory
Infections



Neonatal
Testing



Gastrointestinal
Infections



Bloodstream
Infections



Women's
Health



Healthcare-Associated
Infections



ARIES



LIAISON[®] MDX



VERIGENE[™]

Visit luminexcorp.com and molecular.diasorin.com to learn more.

IVD - For In Vitro Diagnostic Use. Products are region specific and may not be approved in some countries/regions. Please contact Luminex at support@luminexcorp.com to obtain the appropriate product information for your country of residence.

©2022 Luminex Corporation. A DiaSorin Company. All rights reserved. ARIES, VERIGENE, MAGPIX, xTAG, NxTAG, Luminex, and xMAP are trademarks of Luminex Corporation, registered in the US and other countries. Luminex 200 is a trademark of Luminex Corporation. Liaison and Simplexa are trademarks of DiaSorin Molecular, registered in the US and other countries.

KRONUS is Pleased to Offer the First Commercially Available 3-Screen Islet Cell Autoantibody ELISA Kit*

This new 3-Screen ELISA allows for the simultaneous detection and measurement of autoantibodies to GAD and/or IA-2 and/or ZnT8 in a simple, highly-robust, user-friendly assay format.

To obtain additional information regarding this new Product, together with KRONUS' complete offering of immunoassay test kits, please call us toll-free at **800 4 KRONUS** or email us at **kronus@kronus.com**.



ISO 13485 : 2016 QMS Certified

Your Source for Sensitive
Autoimmune Diagnostics

800 4 KRONUS
www.kronus.com

* For Research Use Only. Not For Use in Diagnostic Procedures.